



## Cole-Parmer® Gyrotory Rocking Shakers

Cole-Parmer®

These slow speed gyrotory motion shakers delicately swirl solutions, which is best for low foaming requirements such as delicate cell culturing, DNA extractions, and histological staining procedures. Platform lip and non slip mat hold vessels securely. Choose smaller capacity for personal use or tight spaces. The larger capacity shakers suit larger laboratories and high throughput applications.

- Adjust tilt angle for optimal mixing
- Built in timer and digital speed display provide reproducible results
- Speed range: 5 – 70RPM
- Shaking angle: 3 – 12°

Code	Alt Ref	Model	Platform Dims, mm	Max Load, kg	Dims, w x d x h, mm	Price
<b>MIX2090</b>	51900-26	SSM3	235 x 235	2.9	240 x 300 x 150	<b>£710.00</b>
<b>MIX2072</b>	51900-28	SSL3	355 x 355	10	360 x 420 x 170	<b>£916.00</b>

### Accessories

Code	Alt Ref	Description	Price
<b>MIX2098</b>	51700-94	Tier system (2 platforms + 8 bars)	<b>£214.00</b>



## Grant-Bio PS-3D 3D Rocker

grant<sup>o</sup>

The PS-3D is small but highly functional, with a low profile and small footprint, fitting neatly into the smallest workspace. It provides a smooth 3D orbital motion for mixing in tubes, culture flasks, dishes and boxes, the gentle movement is ideal for staining and destaining of gels without destroying the gel edges, or for keeping fragile tissue intact during incubations. The 3D motion provides very effective sample coverage, for efficient mixing with a reduced quantity of solution. The reliable motor produces regulated and reproducible rotating throughout the speed range and, combined with the sturdy construction, will guarantee years of consistent operation.

- Simple speed adjustment
- Robust drive mechanism provides enhanced reliability – capable of 7 days continuous operation
- Smooth, reliable, extremely quiet motor; low power consumption
- Smooth, non slip mat supplied as standard – prevents vessels from slipping
- Low voltage power supply provides for safe cold room operation and low energy consumption

Code	Alt Ref	Speed, oscil/min	Tilt Angle, °	Max Load, kg	Platform Dims, mm	Dims, w x d x h, mm	Weight, kg	Price
<b>MIX4130</b>	PS-3D	5 - 60	7 (fixed)	1	210 x 210	235 x 235 x 140	1.2	<b>£338.00</b>



### Accessories

Code	Alt Ref	Description	Price
<b>MIX4104</b>	PDM	Dimpled mat to prevent tubes from rolling	<b>£41.00</b>