

## Corning® LSE™ Platform Rockers



A compact yet high capacity rocker that is ideal for Western blotting and general mixing applications. The adjustable rocking speed and fixed 7 degree tilt makes it optimal for staining fragile gels, as well as the more intensive action necessary for washing blots. The speed can be adjusted from 2 to 60RPM and the timer can be set for up to 2 hours or for continuous operation. The unit is safe to use in temperature controlled environments up to 70°C and is available in single or high capacity double platform configurations.

- Low profile and compact size
- 5kg load capacity
- Non slip rubber pads to prevent sample containers from moving during operation
- Double platform model available for increased capacity
- Timer: 0 to 2 hours or continuous
- Speed: 5 to 120 tilts/minute
- Platform size: 300 x 300mm
- Dims, w x d x h: 325 x 360 x 145mm



Code	Description	Price
<b>6704</b>	Single platform	<b>£770.00</b>
<b>6706</b>	Double platform	<b>£843.00</b>

## Thermo Scientific™ Rocker



Rocker with digital display, big platform, accurate speed control and quick angle adjustment for easy sample agitation. Designed for ELISA, protein synthesis, gel staining and destaining, antibody staining, and other time sensitive applications.

- Digital display of speed and time
- Precise speed control:  $\pm 1 \leq 50\text{RPM}$ ,  $\pm 2 > 50\text{RPM}$
- Tray size fits 6 x well plates
- Timer range: continuous or 1min to 99hrs 59min
- Includes dimpled rubber mat and 6 x rubber strips



Code	Alt Ref	Speed Range, RPM	Tilt Angle, °	Max Load, kg	Platform Dims, mm	Dims, w x d x h, mm	Price
<b>MIX1304</b>	88882002	1 - 70	1 - 15	4.5	282 x 322 x 19	387 x 487 x 354	<b>£958.00</b>

## Accessories

Code	Alt Ref	Description	Pack	Price
<b>MIX1314</b>	88882101	Rubber mat, 322 x 282mm	Each	<b>£38.00</b>
<b>MIX1316</b>	88882102	Dimpled rubber mat, 322 x 282mm	Each	<b>£38.00</b>
<b>MIX1318</b>	88882103	Rubber strips	6	<b>£8.00</b>
<b>MIX1352</b>	88882128	Replacement rubber tube	Each	<b>£42.00</b>
<b>MIX1354</b>	88882129	Screws for rubber strips	8	<b>£18.00</b>