



Thermo Scientific™ Mini Rotator

Thermo
SCIENTIFIC

Provides smooth, accurate and orbital mixing of various vessels in the lab. Designed for culturing bacteria and fungi, hybridisation assays, and gel staining and destaining.

- Digital display
- Precise speed control: $\pm 1 \leq 100\text{RPM}$, $\pm 1\% > 100\text{RPM}$
- Continuous operation or timer modes
- Maximum load up to 4kg
- Includes universal aluminium platform

Code	Alt Ref	Speed Range, RPM	Orbit, mm	Platform Dims, mm	Dims, w x d x h, mm	Price
MIX1310	88882008	50 - 500	10	284 x 310	387 x 487 x 264	£1214.00

Accessories

Code	Alt Ref	Description	Pack	Price
MIX1320	88882104	Universal attachment	Each	£217.00
MIX1322	88882105	Universal platform, aluminium	Each	£105.00
MIX1324	88882106	Flask clamp, 50mL (max 20 per platform)	Each	£9.00
MIX1326	88882107	Flask clamp, 100/125mL (max 8 per platform)	Each	£10.00
MIX1328	88882108	Flask clamp, 250mL (max 8 per platform)	Each	£13.00
MIX1330	88882109	Flask clamp, 500mL (max 6 per platform)	Each	£18.00
MIX1332	88882110	Flask clamp, 1000mL (max 4 per platform)	Each	£24.00
MIX1334	88882111	Flask clamp, 2000mL (max 2 per platform)	Each	£31.00
MIX1336	88882112	Sticky mat, 140 x 140mm	Each	£43.00
MIX1338	88882113	Flat platform with rubber mat	Each	£113.00
MIX1346	88882117	Tube rack for 1.5mL tubes	Each	£65.00
MIX1348	88882118	Tube rack for 15mL tubes	Each	£77.00
MIX1350	88882119	Tube rack for 50mL tubes	Each	£77.00
MIX1356	88882130	Screw, GB/T819.1, 4 x 8	4	£6.00
MIX1358	88882131	Screw, GB/T819.1, 5 x 10	4	£6.00
MIX1362	88882133	Screw, GB/T9074.4, 5 x 10	6	£6.00



Thermo Scientific™ Waving Rotator

Thermo
SCIENTIFIC

360° orbital waving rotator for efficient sample mixing with accurate speed control and easy angle adjustment. Designed for ELISA, protein synthesis, gel staining and destaining gels, antibody staining, and other time sensitive applications.

- Digital display of speed and time
- Precise speed control: $\pm 1 \leq 30\text{RPM}$, $\pm 2 > 30\text{RPM}$
- Max load up to 9.8kg
- Tray size fits 6 x well-plates
- Timer range: 1min to 99hrs 59min
- Includes dimpled rubber mat and 6 x rubber strips

Code	Alt Ref	Speed Range, RPM	Tilt Angle, °	Platform Dims, mm	Dims, w x d x h, mm	Price
MIX1306	88882004	5 - 120	1 - 12	276 x 316 x 19	387 x 487 x 354	£1202.00

Accessories

Code	Alt Ref	Description	Pack	Price
MIX1314	88882101	Rubber mat, 322 x 282mm	Each	£38.00
MIX1316	88882102	Dimpled rubber mat, 322 x 282mm	Each	£38.00
MIX1318	88882103	Rubber strips	6	£8.00
MIX1352	88882128	Replacement rubber tube	Each	£42.00
MIX1354	88882129	Screws for rubber strips	8	£18.00