

## Eppendorf Innova® S44i Biological Shakers Continued

eppendorf

## Platforms

| Code    | Alt Ref    | Description        | Capacity                     | Price    |
|---------|------------|--------------------|------------------------------|----------|
| SHA3072 | S44I040001 | Universal platform | —                            | £1122.00 |
| SHA3074 | S44I040002 | Dedicated platform | 86 x 125mL Erlenmeyer flasks | £1465.00 |
| SHA3076 | S44I040003 | Dedicated platform | 53 x 250mL Erlenmeyer flasks | £1354.00 |
| SHA3078 | S44I040004 | Dedicated platform | 39 x 500mL Erlenmeyer flasks | £1271.00 |
| SHA3080 | S44I040005 | Dedicated platform | 23 x 1L Erlenmeyer flasks    | £1203.00 |
| SHA3082 | S44I040006 | Dedicated platform | 15 x 2L Erlenmeyer flasks    | £1136.00 |
| SHA3084 | S44I040007 | Dedicated platform | 8 x 2.8L Fernback flasks     | £1203.00 |
| SHA3086 | S44I040008 | Dedicated platform | 8 x 4L Erlenmeyer flasks     | £1136.00 |



## Eppendorf Excella™ Incubated/Refrigerated Shakers

eppendorf

4 models to suit a wide variety of culturing applications. E24/E24R are benchtop models and accommodate flasks up to 2.8L. E25/E25R are floor standing and accommodate flasks up to 6L. Available with a large selection of accessories to accommodate a wide range of flasks, test tubes and general labware. Like most NBS shakers, Excella™ platforms, clamps and accessories are interchangeable with earlier NBS models.

- Temperature range from 7°C above ambient to 60°C
- Refrigerated versions (R) offer a temperature range of 15°C below ambient to 60°C (minimum set point is 4°C)
- Shaker speed is from 50 to 400RPM
- Platforms are sold separately
- **Note:** see page 680 for platforms – Innova® 2100 platforms are the same ones used in the E24/E24R, and Innova® 2300 platforms are the same ones used in the E25/E25R

| Code    | Alt Ref    | Model | Description    | Orbit, mm | Platform Dims, mm | Dims, w x d x h, mm | Price     |
|---------|------------|-------|----------------|-----------|-------------------|---------------------|-----------|
| SHA4306 | M1352-0002 | E24   | Benchtop       | 19        | 460 x 460         | 560 x 762 x 610     | £6883.00  |
| SHA4307 | M1352-0006 | E24R  | Benchtop       | 19        | 460 x 460         | 560 x 762 x 610     | £9382.00  |
| SHA4308 | M1353-0002 | E25   | Floor standing | 25        | 460 x 760         | 1270 x 770 x 1059   | £11777.00 |
| SHA4310 | M1353-0006 | E25R  | Floor standing | 25        | 460 x 760         | 1270 x 770 x 1059   | £14146.00 |



## Eppendorf Innova® Incubated/Refrigerated Shakers

eppendorf

Heavy duty triple eccentric drive providing uniform motion and superb reliability. All units have audible and visual alarms, speed 25-500RPM, timer (0-99.9hrs) with continuous operation or automatic shut off at the end of the study.

- Temperature range is from 5°C above ambient to 80°C
- Refrigerated versions (R) offers a temperature range of 20°C below ambient to 80°C (minimum set point is 4°C)
- 40/40R models have been designed for benchtop use with a maximum flask size of 3L
- 43/43R models are floorstanding and can accommodate 6L flasks
- Platforms are sold separately
- **Note:** see page 680 for platforms – Innova® 2100 platforms are the same ones used in the 40/40R and Innova® 2300 platforms are the same ones used in the 43/43R

| Code    | Alt Ref    | Model | Orbit, mm | Dims, w x d x h, mm | Platform Dims, mm | Price     |
|---------|------------|-------|-----------|---------------------|-------------------|-----------|
| SHA4131 | M1299-0082 | 40    | 19        | 559 x 762 x 610     | 460 x 460         | £8927.00  |
| SHA4130 | M1299-0092 | 40    | 25        | 559 x 762 x 610     | 460 x 460         | £8927.00  |
| SHA4133 | M1299-0086 | 40R   | 19        | 559 x 762 x 610     | 460 x 460         | £10472.00 |
| SHA4132 | M1299-0096 | 40R   | 25        | 559 x 762 x 610     | 460 x 460         | £10472.00 |
| SHA4101 | M1320-0012 | 43    | 51        | 1156 x 813 x 1031   | 760 x 460         | £14956.00 |
| SHA4100 | M1320-0002 | 43    | 25        | 1156 x 813 x 1031   | 760 x 460         | £14956.00 |
| SHA4111 | M1320-0016 | 43R   | 51        | 1156 x 813 x 1031   | 760 x 460         | £16955.00 |
| SHA4110 | M1320-0006 | 43R   | 25        | 1156 x 813 x 1031   | 760 x 460         | £16955.00 |