

Grant JB Academy Unstirred Water Baths



A basic, simple to use, high quality water bath range, consisting of three great value models. A polycarbonate, perforated base tray is included as standard (lid not included), which promotes heat convection, optimal temperature uniformity, and allows any container to be placed directly in the bath.

- Temperature range: ambient +5 to 95°C ($\pm 0.5^\circ\text{C}$)
- Unique Set and Forget® technology with fast heat up
- Practical front panel lock which disables front panel controls preventing unintentional temperature changes



Code	Alt Ref	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x l, mm	Drain Outlet	Price
BAT3310	JBA5	5	335 x 215 x 205	281 x 131	No	£389.00
BAT3312	JBA12	12	360 x 390 x 227	306 x 281	No	£448.00
BAT3314	JBA18	18	335 x 590 x 277	281 x 485	No	£520.00

Grant Lids for Unstirred Water Baths



Use lids to direct condensation away from immersed vessels, avoid contamination, reduce evaporation and save energy. Recommended when using water baths above 60°C. Flat lids have ring sets of variable hole diameter to accommodate tall vessels.



Code	Description	For use with	Price
BAT3031	Gabled polycarbonate lid	SAP2	£58.00
BAT3036	Gabled polycarbonate lid	JBN5, JBA5, SAP2S, SAP5, SBBAQUA5+	£73.00
BAT3037	Gabled polycarbonate lid	JBN12, JBA12, SAP12, SBBAQUA12+	£89.00
BAT3039	Gabled polycarbonate lid	JBN18, JBN26, JBA18, SAP18, SAP26, SBBAQUA26+	£98.00
BAT3120	Gabled stainless steel lid	JBN5, JBA5, SAP2S, SAP5, SBBAQUA5+	£138.00
BAT3122	Gabled stainless steel lid	JBN12, JBA12, SAP12, SBBAQUA12+	£151.00
BAT3124	Gabled stainless steel lid	JBN18, JBN26, JBA18, SAP18, SAP26, SBBAQUA18+, SBBAQUA26+	£173.00
BAT3126	Gabled stainless steel lid	SAP34	£199.00
BAT3130	Flat stainless steel lid with 2 ring sets	JBN5, JBA5, SAP5, SBBAQUA5+	£167.00
BAT3132	Flat stainless steel lid with 4 ring sets	JBN12, JBA12, SAP12, SBBAQUA12+	£256.00
BAT3134	Flat stainless steel lid with 6 ring sets	JBN18, JBN26, JBA18, SAP18, SAP26, SBBAQUA18+, SBBAQUA26+	£362.00
BAT3136	Flat stainless steel lid with 8 ring sets	SAP34	£418.00