

Genlab Large Capacity Ovens Continued

Horizontal Style with Twin Doors



Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Int Dims, w x d x h, mm	Ext Dims, w x d x h, mm	Price
OVE2008	LCO/250H/TDIG	220	3	4	920 x 440 x 540	1090 x 570 x 940	£3447.00
OVE2010	LCO/350H/TDIG	320	3	4	920 x 640 x 540	1090 x 770 x 940	£3870.00
OVE2012	LCO/425H/TDIG	450	4	5	920 x 640 x 760	1440 x 810 x 1120	£4747.00
OVE2014	LCO/500H/TDIG	535	5	7	920 x 640 x 910	1440 x 810 x 1270	£5135.00
OVE2016	LCO/650H/TDIG	645	5	7	920 x 770 x 910	1600 x 950 x 1270	£7302.00
OVE2018	LCO/750H/TDIG	750	5	8	920 x 770 x 1060	1600 x 950 x 1420	£7541.00
OVE2020	LCO/850H/TDIG	865	5	9	920 x 770 x 1220	1600 x 950 x 1580	£7791.00
OVE2022	LCO/1000H/TDIG	1005	5	11	920 x 770 x 1420	1600 x 950 x 1780	£8710.00
OVE2024	LCO/1250H/TDIG	1280	5	11	1220 x 740 x 1420	1900 x 920 x 1780	£9748.00

Genlab Dual Purpose Ovens/Incubators



A range of dual purpose ovens/incubators, which combine the advantage of two temperature ranges in the same unit. Using bespoke firmware designed specifically for this chamber size, and a calibration class temperature sensor, we can offer accurate control up to 250°C without temperature overshoot. The same chamber can also maintain a steady 30°C constant temperature by simply adjusting the set point. Available in a variety of sizes and designed with both side and base mounted displays to suit your exact needs.

- Temp range: ambient +8 to 250° C
 - Resolution: 0.1°C
 - On screen historical trending (48 hours)
 - Manual or automatic overheat reset
 - High accuracy Pt100A duplex sensors <0.6°C
- Low chamber uniformity
 - Energy efficient with integral energy meter
 - Set point lock feature



Base Mount Design

Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
INC2242	MINIS/6/TDIG	6	1	1	350 x 330 x 460	230 x 190 x 150	£940.00
INC2244	MINIS/18/TDIG	18	1	1	380 x 410 x 550	260 x 270 x 260	£1023.00
INC2246	MINIS/30/TDIG	30	2	2	480 x 490 x 510	360 x 350 x 240	£1096.00
INC2248	MINIS/40/TDIG	40	2	3	480 x 490 x 540	360 x 350 x 320	£1208.00
INC2250	MINIS/50/TDIG	53	2	3	630 x 660 x 470	490 x 330 x 330	£1279.00
INC2252	MINIS/75/TDIG	73	2	3	630 x 660 x 600	490 x 450 x 330	£1480.00