

Astell Compact Autoclaves Continued



Options and Accessories

Code	Alt Ref	Description	Price
AUT4310	AAR130	Integral data printer	£1329.00
AUT4182	AAR014	Load sensed process timing	£555.00
AUT4418	AAR122	Ethernet interface	£323.00
AUT4018	AAN346	Discard container, each model holds 33L (1), 43L (2), 63L (3)	£478.00
AUT4320	AAN025	5 position shelf rack (plus 2 shelves) for AUT4300/AUT4400	£671.00
AUT4322	AAN022	Spare shelf for AUT4300/AUT4400	£285.00
AUT4324	AAN530	5 position shelf rack (plus 3 shelves) for AUT4220/AUT4402	£735.00
AUT4326	AAN532	Spare shelf for AUT4220/AUT4402	£310.00
AUT4328	AAN040	5 position shelf rack (plus 3 shelves) for AUT4304/AUT4404	£762.00
AUT4330	AAN043	Spare shelf for AUT4304/AUT4404	£362.00
AUT4016	AAN308	Spare condensate bottle for Classic models	£207.00
AUT4253	AAN007	Air assisted cooling	£336.00

Astell Compact Autoclaves with Closed Door Drying



Designed to offer convenience and flexibility in the sterilisation of wrapped and unwrapped instruments. This range of 'vacuum-less drying' steam sterilisers include Astell's 'Ecofill' water management system as standard.

- Four programmed cycles suitable for solid metal instruments (unwrapped), solid metal instruments (wrapped), glassware (unwrapped), and hand pieces in surgical pouches
- Closed door drying - following sterilisation, filtered air is circulated through the chamber by an efficient air pump rapidly drying each load
- Included furniture: chamber shelf rack, two wire shelves and one pouch rack
- Multi-programmable touch screen controller
- Front loading, cylindrical chambers



Code	Alt Ref	Capacity, L	Chamber Dims, Ø x d, mm	Overall Dims, w x d x h, mm	Price
AUT2068	AMB220D	33	350 x 360	520 x 805 x 600	£11543.00
AUT2078	AMB230D	43	350 x 460	520 x 805 x 600	£13036.00
AUT2088	AMB240D	63	350 x 660	520 x 1005 x 600	£14905.00

Sustainable laboratory solutions!

Visit the SLS website to read online or request your copy:

scientificlabs.co.uk/sustainability (UK & Northern Ireland)

scientificlabs.ie/sustainability (Republic of Ireland)

scientificlabs.com/sustainability (International)

