

Nickel Electro Clifton® NE1 Unstirred Water Baths



Designed for general use and are economically priced without compromising quality. With a high quality corrosion resistant stainless steel tank, all painted surfaces have an antibacterial finish with good chemical resistance. Heater and temperature sensors are concealed under the tank.

- Analogue dial for setting temperature, graduated in 5°C divisions
- Sensitivity: $\pm 0.25^{\circ}\text{C}$
- Uniformity: $\pm 0.1^{\circ}\text{C}$
- Temperature range: ambient +5 to 100°C (not designed for continuous boiling)
- Heater with safety cut-out
- Lids recommended above 60°C (sold separately)
- Includes shallow height stainless steel perforated false base
- Drain outlet on 22 and 28L baths, option of accessory gated tap 4



Code	Alt Ref	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Power, W	Drain Outlet	Price
BAT7000	NE1-4	4	332 x 185 x 290	300 x 150 x 150	500	No	£349.00
BAT7002	NE1-8	8	332 x 270 x 290	300 x 240 x 150	1000	No	£385.00
BAT7005	NE1-14	14	361 x 332 x 290	325 x 300 x 150	1000	No	£409.00
BAT7010	NE1-22	22	538 x 332 x 290	500 x 300 x 150	1500	Yes	£487.00
BAT7012	NE1-28	28	538 x 332 x 290	500 x 300 x 200	1500	Yes	£531.00

Nickel Electro Clifton® NE1D Digital Unstirred Water Baths



Clifton digital baths are designed for general use and are economically priced without compromising quality. Simple digital controller for setting temperature with 1°C divisions. Designed for immersion of flasks, bottles and racks in a highly stable temperature environment.

- Display/resolution: digital LED/ 1°C
- Temperature range: ambient +5 to 99°C
- Temperature setting: push button
- Stability/uniformity: $\pm 0.2^{\circ}\text{C}/\pm 0.1^{\circ}\text{C}$
- Temperature memory retained even after power off
- Complete with perforated shelf for even heat distribution
- Lids recommended above 60°C (supplied separately)



Code	Alt Ref	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Power, W	Price
BAT2838	NE1D-4	4	332 x 185 x 290	300 x 150 x 150	400	£381.00
BAT2839	NE1D-8	8	332 x 270 x 290	300 x 240 x 150	800	£401.00
BAT2840	NE1D-14	14	361 x 332 x 290	325 x 300 x 150	1000	£446.00
BAT2842	NE1D-22	22	538 x 332 x 290	500 x 300 x 150	1500	£539.00

Can't find what you're looking for?

Let us know! Our team are available to help. Contact your local Territory Sales Manager to discuss your requirements.



**SCIENTIFIC
LABORATORY
SUPPLIES**

