



## Nickel Electro Clifton® NE2D Digital PID Controlled Unstirred Water Baths



Provides a stable temperature environment with digital PID temperature control for precision. Controls feature auto alarm settings, illuminated on/off switch, heater and over temperature indications.

- Sensitivity:  $\pm 0.2^{\circ}\text{C}$
- Uniformity:  $\pm 0.1^{\circ}\text{C}$
- Temperature range: ambient +5 to  $99^{\circ}\text{C}$
- Digital PID control with LED temperature display
- Over temperature alarm with heater cut off  $+4^{\circ}\text{C}$  above set point
- Touch pad controls with easy wipe surfaces
- Painted body and controls features antibacterial protection, hygienic coating which actively inhibits bacterial growth
- Lid recommended above  $60^{\circ}\text{C}$  (sold separately)
- 2 and 5L baths are shallow

| Code           | Alt Ref | Capacity, L | Ext Dims, w x d x h, mm | Int Dims, w x d x h, mm | Power, W | Drain Outlet | Price          |
|----------------|---------|-------------|-------------------------|-------------------------|----------|--------------|----------------|
| <b>BAT7038</b> | NE2-2D  | 2           | 332 x 185 x 135         | 300 x 150 x 65          | 200      | No           | <b>£498.00</b> |
| <b>BAT7022</b> | NE2-4D  | 4           | 332 x 185 x 290         | 300 x 150 x 150         | 400      | No           | <b>£532.00</b> |
| <b>BAT7035</b> | NE2-5D  | 5           | 361 x 332 x 135         | 325 x 300 x 65          | 800      | No           | <b>£509.00</b> |
| <b>BAT7024</b> | NE2-8D  | 8           | 332 x 270 x 290         | 300 x 240 x 150         | 800      | No           | <b>£504.00</b> |
| <b>BAT7036</b> | NE2-10D | 10          | 332 x 270 x 290         | 300 x 240 x 200         | 800      | No           | <b>£650.00</b> |
| <b>BAT7027</b> | NE2-14D | 14          | 361 x 332 x 290         | 325 x 300 x 150         | 1000     | No           | <b>£580.00</b> |
| <b>BAT7032</b> | NE2-22D | 22          | 538 x 332 x 290         | 500 x 300 x 150         | 1500     | Yes          | <b>£618.00</b> |
| <b>BAT7034</b> | NE2-28D | 28          | 538 x 332 x 290         | 500 x 300 x 200         | 1500     | Yes          | <b>£745.00</b> |



## Nickel Electro Clifton® Lids for Unstirred Water Baths



Lids help to reduce water evaporation, assists in keeping samples and water free from contamination and saves energy. Clear polycarbonate gabled lids can be hinged back out of the way to access the water bath. Stainless steel gabled lids have a large insulated lifting handle.

| Code             | Alt Ref  | Description                     | For use with                                                                                             | Price          |
|------------------|----------|---------------------------------|----------------------------------------------------------------------------------------------------------|----------------|
| <b>BAT7070</b>   | GL1-4    | Hinged polycarbonate gabled lid | NE2-4D, NE1-4, NE1D-4                                                                                    | <b>£73.00</b>  |
| <b>BAT7071</b>   | GL1-8    | Hinged polycarbonate gabled lid | NE2-10D, NE2-8D, NE1-8, NE5-10D                                                                          | <b>£73.00</b>  |
| <b>BAT7072</b>   | GL1-14   | Hinged polycarbonate gabled lid | NE1D-14, NE1B-18W, NE1-14, NE2-14D, NE1B-14                                                              | <b>£73.00</b>  |
| <b>BAT7160</b>   | GL1-22   | Hinged polycarbonate gabled lid | NE2-28D, NE5-28, NE1-22, NE2-22D, NE1-28, NE1B-28, NE5-28D, NE1D-22                                      | <b>£73.00</b>  |
| <b>BAT7090</b>   | SL1-4    | Stainless steel gabled lid      | NE1-2/2/2, NE2D-4/4, NE1D-4/4/4, NE2-2D, NE1D-2/2/2, NE3-4DT, NE1-4/4/4, NE2-4D, NE2D-4/8, NE1-4, NE1D-4 | <b>£134.00</b> |
| <b>BAT7163</b>   | SL1-8    | Stainless steel gabled lid      | NE1-8, NE5-10D, NE2-8D, NE2-10D, NE2D-8/22, NE2D-4/8, NE1D-8, NE3-8DT                                    | <b>£143.00</b> |
| <b>BAT7092</b>   | SL1-14   | Stainless steel gabled lid      | NE2-5D, NE3-14DT, NE2-14D, NE1-14, NE1D-14                                                               | <b>£143.00</b> |
| <b>BAT7094</b>   | SL1-22   | Stainless steel gabled lid      | NE5-28D, NE2-28D, NE1-9, NE1D-22, NE1-28, NE2D-8/22, NE2-9D, NE3-28DT, NE5-28, NE2-22D, NE1-22           | <b>£159.00</b> |
| <b>BAT7162</b>   | SL1-22H  | Stainless steel gabled lid      | NE5-28D, NE1D-22, NE2-28D, NE1-28, NE3-28DT, NE5-28, NE2-22D, NE1-22                                     | <b>£273.00</b> |
| <b>EDGL14AF</b>  | GL1-4AF  | Clear polycarbonate gabled lid  | All 4L unstirred baths                                                                                   | <b>£73.00</b>  |
| <b>EDGL18AF</b>  | GL1-8AF  | Clear polycarbonate gabled lid  | All 8L unstirred baths                                                                                   | <b>£73.00</b>  |
| <b>EDGL114AF</b> | GL1-14AF | Clear polycarbonate gabled lid  | All 14L unstirred baths                                                                                  | <b>£73.00</b>  |
| <b>EDGL122AF</b> | GL1-22AF | Clear polycarbonate gabled lid  | All 22L unstirred baths                                                                                  | <b>£73.00</b>  |