

Eppendorf HeatSealer Microplate Sealer

eppendorf

The Eppendorf HeatSealer Family offers you reliable and effective sealing of a variety of plates. You can choose from a more compact (S100) and a more flexible (S200) instrument according to your sealing needs.

- Hermetic sealing of multi-well plates
- Minimises evaporation in PCR, reducing cross contamination
- Seals 96- and 384-well plates
- Integrated thermostat prevents overheating
- Improved mechanics for easier sealing
- PCR plate adapter included



| Code | Model | Temperature, °C | Weight, kg | Dims, w x d x h, mm | Power, W | Price |
|-------------|-------|-----------------|------------|---------------------|----------|----------|
| E5391000036 | S100 | 170 | 4 | 220 x 250 x 230 | 320 | £2302.00 |
| E5392000030 | S200 | 125 to 200 | 7.2 | 220 x 321 x 425 | 350 | £3304.00 |

Accessories

| Code | Description | For use with | Price |
|-------------|----------------------------|--------------|---------|
| E5391070018 | Low profile plate adapter | S100 | £294.00 |
| E5392070020 | Low profile plate adapter | S200 | £294.00 |
| E5391070034 | High profile plate adapter | S100 | £294.00 |
| E5392070011 | High profile plate adapter | S200 | £294.00 |
| E5391070026 | 96-well PCR plate adapter | S100 | £294.00 |
| E5392070038 | 96-well PCR plate adapter | S200 | £294.00 |


Eppendorf Heat Sealing Foil/Film

eppendorf

- Hermetic sealing of multi-well plates, especially recommended for low reaction volumes
- Best protection against evaporation during PCR
- Maximum adhesiveness at high temperatures during PCR

Supplied in packs of 100

| Code | Description | Price |
|-------------|-------------------|---------|
| E0030127854 | Heat sealing foil | £114.00 |
| E0030127838 | Heat sealing film | £162.00 |



BRAND 0.2mL PCR Tubes



- Free of DNase, DNA, RNase, and endotoxins for reliable and contaminant free results
- Thin walls for rapid heat transfer
- Tightly closing covers protect samples against evaporation

Supplied in packs of 1000

| Code | Alt Ref | Description | Price |
|---------|---------|-------------|--------|
| MOL2300 | 781300 | Domed cap | £73.00 |
| MOL2302 | 781305 | Flat cap | £60.00 |

