


SLS Select Petri Dishes



- High quality, disposable Petri dishes
- Aseptically manufactured in clear polystyrene to the highest standards
- Compatible with most automatic systems



Code	Description	Width, mm	Pack	Price
PET3000	Circular, no vent	55	1620	£214.04
PET3003	Circular contact plate, triple vent, grip lid	65	720	£187.16
SLS2000	Circular, single vent	90	500	£87.00
SLS2002	Circular, triple vent	90	500	£87.00
PET3004	Circular, no vent	90	825	£133.39
PET3005	Circular, 3 compartment, triple vent	90	552	£156.00
P903S720	Circular, low style, triple vent	90	720	£175.00
PET3008	Square, four vents	120	252	£169.58
PET3007	Circular, triple vents	140	120	£125.00



Thermo Scientific™ Sterilin™ Standard Petri Dishes



- Shallow dishes maximise incubator space
- Manufactured from non-cytotoxic virgin polystyrene
- For in vitro use only
- 90mm dishes are ideal for use in automatic plate pourers
- 50mm deep form dish is over 20mm deep and designed for use with liquid media to aid handling
- 140mm dish is produced and tested in accordance with the EN ISO 24998 Petri dish standard which includes stringent dimensional controls

Available in the following formats:

- Non vented - most suitable for anaerobic and long-term work
- Single vent - limits gaseous exchange, minimizes evaporation and dehydration; ideally suited for long-term work
- Triple vent - aids gaseous exchange; ideally suited for short-term work

Code	Alt Ref	Description	Dims, Ø x h, mm	Sterility	Inner Pack	Pack	Price
PET2038	101RT/C	Non vented	90 x 15.7	Aseptically manufactured	20	500	£136.00
PET2040	101RT/IRR	Non vented	90 x 15.7	Sterile	20	500	£352.00
PET2020	122-17	Single vent	50 x 14	Aseptically manufactured	10	700	£365.00
PET2022	124-17	Single vent	50 x 20.3	Aseptically manufactured	20	500	£375.00
PET2026	123-17	Single vent	60 x 15	Aseptically manufactured	10	540	£342.00
PET2000	101R20	Single vent	90 x 15.9	Aseptically manufactured	20	500	£136.00
PET2001	101/IRR	Single vent	90 x 15.9	Sterile	20	500	£246.00
PET2018	121V	Triple vent	35 x 11	Aseptically manufactured	10	800	£592.00
PET2036	VR5B	Triple vent	55 x 12	Aseptically manufactured	15	1620	£596.00
PET2002	101VR20	Triple vent	90 x 16.1	Aseptically manufactured	20	500	£136.00
PET2003	101V/IRR	Triple vent	90 x 16.1	Sterile	20	500	£246.00
PET2032	501V	Triple vent	140 x 21.1	Aseptically manufactured	10	80	£174.00