

SLS Select Serological Pipettes



- Manufactured in compliance with the current versions of the ISO 13485 and FDA CFR 21 Part 820, cGMP
- Made from USP Class VI polystyrene resin
- Crisp, dark ascending and descending as well as negative graduations
- Colour-coded rings for easy identification
- Plugged
- Non pyrogenic
- RNase/DNase free
- Sterile



Code	Volume, mL	Packaging	Inner Pack	Pack	Price
SLS4054	1	Individually wrapped	100	1000	£172.00
SLS4056	1	Bulk	25	1000	£135.00
SLS4058	2	Individually wrapped	100	1000	£184.00
SLS4060	2	Bulk	25	1000	£134.00
SLS4064	5	Individually wrapped	50	200	£47.50
SLS4062	5	Bulk	25	500	£125.00
SLS4066	10	Individually wrapped	50	200	£49.90
SLS4068	10	Bulk	25	500	£161.00
SLS4072	25	Individually wrapped	50	200	£112.00
SLS4070	25	Bulk	20	200	£187.00
SLS4074	50	Individually wrapped	20	100	£144.00

Thermo Scientific™ Sterilin™ Shortie Serological Pipettes



- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level <0.01EU/mL
- Batch-specific certificate supplied in each carton containing lot number for complete traceability
- Accuracy of $\pm 1\%$ for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting, with negative graduations for extra capacity
- Colour-coded for ease of identification
- Validated non haemolytic
- Supplied with suction adapter, unless otherwise stated

Supplied in packs of 200



Code	Alt Ref	Capacity, mL	Negative Graduations, mL	Price
PIP4016	42505	5	-3.0	£260.00
PIP4028	42510	10	-10.0	£386.00