

Thermo Scientific™ Sterilin™ Serological Pipettes

Thermo
SCIENTIFIC

- Gamma irradiated to ensure sterility
- Certified non-pyrogenic with an endotoxin level <0.01EU/mL
- Batch-specific certificate supplied in each carton containing lot number for traceability
- Accuracy of $\pm 1\%$ for nominal volume
- Clear black printed graduations for maximum clarity
- Descending and ascending graduations to aid pipetting, with negative graduations for extra capacity
- Colour coded for ease of identification
- Validated non haemolytic



Code	Alt Ref	Capacity, mL	Description	Inner Pack	Pack	Price
PIP4001	40101K	1	Plugged, individually wrapped	100	1000	£306.00
PIP4003	40301K	1	Plugged, bulk	50	1000	£289.00
PIP4005	41301K	1	Open ended, bulk	50	1000	£298.00
PIP4007	40102K	2	Plugged, individually wrapped	100	1000	£355.00
PIP4010	40105	5	Individually wrapped	1	200	£114.00
PIP4015	47105	5	Plugged, individually wrapped	50	200	£128.00
PIP4017	47105N	5	Narrow, plugged, bulk	50	200	£135.00
PIP4020	47110	10	Plugged, individually wrapped	50	500	£319.00
PIP4024	47310	10	Plugged, bulk	25	500	£317.00
PIP4029	47110N	10	Narrow, plugged, individually wrapped	50	500	£319.00
PIP4030	40125	25	Plugged, individually wrapped	25	200	£286.00
PIP4032	47150	50	Plugged, individually wrapped	1	50	£196.00

Thermo Scientific™ Nunc™ Serological Pipettes

Thermo
SCIENTIFIC

- Stringent quality specifications and convenient features make these serological pipettes a great choice for daily lab work
- Bold, clear bi-directional graduations
- Colour coded pipettes and packaging for convenient volume selection
- Individually packed pipettes in either paper/plastic or plastic peel for high sterility assurance
- Reduce packaging waste with bulk packed options
- Certified calibration of accuracy and precise graduations ensures dispensing accuracy $\pm 2\%$ for 1 and 2mL and $\pm 1\%$ for 5, 10, 25 and 50mL at full volume



Code	Capacity, mL	Packaging	Inner Pack	Pack	Price
170353N	1	Paper/plastic peel	1	1000	£158.00
170364N	1	Plastic film	1	1000	£144.00
170371N	1	Bulk	50	1000	£130.00
170354N	2	Paper/plastic peel	1	500	£88.00
170365N	2	Plastic film	1	500	£73.00
170372N	2	Bulk	50	500	£64.00
170355N	5	Paper/plastic peel	1	200	£48.00
170366N	5	Plastic film	1	200	£45.00
170373N	5	Bulk	50	200	£110.00
170356N	10	Paper/plastic peel	1	200	£62.00
170367N	10	Plastic film	1	200	£57.00
170374N	10	Bulk	50	200	£131.00
170357N	25	Paper/plastic peel	1	200	£102.00
170368N	25	Plastic film	1	200	£100.00
170375N	25	Bulk	25	200	£100.00
170358N	50	Paper/plastic peel	1	100	£151.00
170369N	50	Plastic film	1	100	£141.00
170376N	50	Bulk	25	100	£143.00