



Heathrow Scientific Pipette Racks



- Angled four shelved compartments keep your 0.1 to 50mL pipettes or cans within easy reach while preventing damage and clutter
- Heavy duty acrylic version provides crystal-clear viewing from any angle
- ABS versions have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- Magnetised ABS version has strong magnets on both sides and mounts securely to most metal surfaces
- Dims, w x d x h: 88 x 291 x 401mm

Code	Alt Ref	Material	Colour	Price
PIP5510	HS20611	Acrylic	Clear	£154.00
PIP5512	HS20612	ABS	Blue	£97.00
PIP5514	HS20613	ABS	Red	£97.00
PIP5518	HS20612M	ABS magnetised	Blue	£97.00



Heathrow Scientific Universal Pipette Carousels



- Holds up to 8 single channel pipettes, 4 multi channel pipettes, or a combination of both
- Inserts on the upper holders can be added or removed to accommodate different pipette widths
- Rotating head and stable weighted base for one handed removal of pipettes
- Ideal for labs with multiple brands of pipettes
- Conveniently rotate stand head to select or store pipettes
- Quick grip handle to move pipette stand (loaded or unloaded)
- Can be cleaned with soap/water or common laboratory disinfectants (not autoclavable)
- Made of ABS
- Dims, base Ø x top Ø x h: 158 x 229 x 408mm

Code	Alt Ref	Colour	Price
PIP3920	HS120480	Blue/green	£122.00
PIP3922	HS120481	Purple/clear	£127.00
PIP3924	HS120499	Clear/rainbow	£127.00



Heathrow Scientific Pipette Tip Bins



- Functional acrylic storage containers for tips, tubes and syringes
- 3mm thick clear acrylic
- Flip-top covers and generous openings allow for easy filling and dispensing
- Non-skid rubber feet

Code	Alt Ref	Description	Dims, l x w x d, mm	Price
ST08620	HS23401	1 compartment	184 x 194 x 232	£82.00
PIP9400	HS23402	2 compartment	150 x 177 x 368	£158.00
PIP9402	HS23403	3 compartment	217 x 170 x 217	£180.00
PIP9404	HS23404	4 compartment	287 x 170 x 217	£145.00