

Thermo Scientific™ Finnpiquette™ F1 Pipettes

Thermo
SCIENTIFIC

A lightweight design that maximises user comfort and reduces strain. Available as 14 single channel models covering 0.2µL to 10mL and 8, 12 and 16 multichannel versions covering 1 to 300µL.

- One-handed volume changes for convenience
- Optional handle strips for your F1 pipette also enable a firmer grip
- 120° adjustable finger rest – for both right- and left-handed operator comfort
- Super blow-out ensures efficient delivery of micro-volumes and prevents capillary action in 50µL models and below
- Secure volume lock

GLP kit 1 includes:

- Single channel 1-10µL, 10-100µL and 100-1000µL pipettes
- F-Series stand
- 3 racks of Finntip Flex tips

GLP kit 2 includes:

- Single channel 0.2-2µL, 2-20µL, 20-200µL and 100-1000µL pipettes
- F-Series stand
- 4 racks of Finntip Flex tips

GLP kit 3 includes:

- Single channel 10-100µL, 100-1000µL and 1-10mL pipettes
- F-Series stand
- 2 racks of Finntip Flex tips and 1 rack of Finntip tips



GLP kit 4 includes:

- Single channel 2-20µL, 20-200µL and 100-1000µL pipettes
- F-Series stand
- 3 racks of Finntip Flex tips

Code	Alt Ref	Channel	Volume Range, µL	Price
PIP4900	4641010N	Single	0.2 - 2	£371.00
PIP4902	4641020N	Single	0.5 - 5	£371.00
PIP4904	4641030N	Single	1 - 10, micro	£371.00
PIP4906	4641040N	Single	1 - 10	£371.00
PIP4908	4641050N	Single	2 - 20, micro	£371.00
PIP4910	4641060N	Single	2 - 20	£371.00
PIP4926	4641130N	Single	5 - 50, micro	£371.00
PIP4928	4641140N	Single	5 - 50	£371.00
PIP4912	4641070N	Single	10 - 100	£371.00
PIP4914	4641080N	Single	20 - 200	£371.00
PIP4916	4641090N	Single	30 - 300	£371.00
PIP4918	4641100N	Single	100 - 1000	£371.00
PIP4920	4641110N	Single	500 - 5000	£386.00
PIP4922	4641120N	Single	1000 - 10000	£386.00
PIP4956	4661000N	8	1 - 10	£975.00
PIP4958	4661010N	8	5 - 50	£975.00
PIP4960	4661020N	8	10 - 100	£975.00
PIP4962	4661030N	8	30 - 300	£975.00
PIP4964	4661040N	12	1 - 10	£1105.00
PIP4966	4661050N	12	5 - 50	£1105.00
PIP4968	4661060N	12	10 - 100	£1105.00
PIP4970	4661070N	12	30 - 300	£1105.00
PIP4972	4661080N	16	1 - 10	£1182.00
PIP4974	4661090N	16	5 - 50	£1182.00
PIP4976	4700850N	GLP kit 1	1 - 10, 10 - 100, 100 - 1000	£959.00
PIP4978	4700860N	GLP kit 2	0.2 - 2, 2 - 20, 20 - 200, 100 - 1000	£1156.00
PIP4980	4700865N	GLP kit 3	10 - 100, 100 - 1000, 1000 - 10000	£977.00
PIP4982	4701060N	GLP kit 4	2 - 20, 20 - 200, 100 - 1000	£959.00