



## Eppendorf Research® plus Move It® Pipettors with Adjustable Tip Spacing

eppendorf

Do you frequently transfer samples between different vessel formats, such as between tubes and plates, and strive for more efficient and safe solutions? The new adjustable tip spacing pipette Move It® accelerates and simplifies synchronous pipetting of a series of samples. A unique design without tubing connections results in increased precision and durability as well as user and sample safety.

- Efficient sample transfers from one vessel type to another.
- Pre-defined settings allow for quick switches back and forth between different formats
- Tip spacing freely selectable between 4.5 and 33mm
- Time saving thanks to rotatable pipette head - 360°
- Tubeless system for reliability and reproducible results - no cable tangle and tubing-porosity, for consistently reliable robustness and precision

Code	Channel	Volume Range, µL	Tip Spacing Range, mm	Price
E3125000150	4	30 - 300	9 - 33	£1200.00
E3125000184	4	120 - 1200	9 - 33	£1200.00
E3125000168	6	30 - 300	9 - 20	£1249.00
E3125000192	6	120 - 1200	9 - 20	£1249.00
E3125000117	8	2 - 20	4.5 - 14	£1288.00
E3125000133	8	10 - 100	4.5 - 14	£1288.00
E3125000176	8	30 - 300	9 - 14	£1300.00
E3125000206	8	120 - 1200	9 - 14	£1300.00
E3125000125	12	2 - 20	4.5 - 9	£1387.00
E3125000141	12	10 - 100	4.5 - 9	£1387.00

## Corning® Lambda™ Plus Pipette Starter Kit

CORNING

All Lambda™ Plus pipettors feature smooth plunger movement and extremely low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button, and a locking mechanism prevents accidental change. Designed to maintain their calibration, even after repeated use and autoclaving. A key is included with each pipettor for easy in-house calibration.

- Factory calibrated to strict ISO 8655 standards
- Fully autoclavable and UV resistant
- Designed for universal fit with all common brands of pipette tips
- Adjustable tip ejector

### Kit includes:

- 4 x Lambda™ Plus pipettes (0.5 - 10µL, 2 - 20µL, 20 - 200µL, 100 - 1000µL)
- 1 x universal carousel rack for 6 single channel pipettes
- 3 x DeckWorks™ pipette tips (10µL, 200µL, 1000µL)
- 1 x Corning® insulated lunch bag free gift

Code	Price
4069	£733.00

