

Sartorius Picus® Electronic Pipettes Continued **SARTORIUS**

Picus® NxT Electronic Pipettes

Offers all the advantages of the Picus® electronic pipettor with added compliance benefits.

- Electronic tip ejection
- RFID capability and password protection
- Repeat blowout

Code	Alt Ref	Channel	Volume Range, µL	Price
PIP3900	LH-745021	Single	0.2 - 10	£736.00
PIP3901	LH-745041	Single	5 - 120	£736.00
PIP3902	LH-745061	Single	10 - 300	£736.00
PIP3903	LH-745081	Single	50 - 1000	£736.00
PIP3904	LH-745101	Single	100 - 5000	£736.00
PIP3905	LH-745111	Single	500 - 10000	£736.00
PIP3906	LH-745321	8	0.2 - 10	£1132.00
PIP3907	LH-745341	8	5 - 120	£1132.00
PIP3908	LH-745361	8	10 - 300	£1132.00
PIP3909	LH-745391	8	50 - 1200	£1132.00
PIP3910	LH-745421	12	0.2 - 10	£1322.00
PIP3911	LH-745441	12	5 - 120	£1322.00
PIP3912	LH-745461	12	10 - 300	£1322.00
PIP3913	LH-745491	12	50 - 1200	£1322.00

Accessories

Code	Alt Ref	Description	Price
PIP7762	730981	Single place charging stand with AC adapter	£226.00
PIP9506	730991	Charging carousel for 4 pipettes	£408.00
PIP7678	725620	Linear stand for up to 12 pipettes	£69.00

Thermo Scientific™ Finnpiquette™ Novus Electronic Pipettes

Thermo
SCIENTIFIC

Achieve improved ergonomics and performance in standard pipetting applications. The user interface simply displays all pipetting actions, and the index-finger trigger pipetting action prevents common Repetitive Strain Injuries (RSI). Simplicity, comfort and functionality.

- Easy-to-read and intuitive menu navigation
- 9 aspirate/dispense speeds
- Save up to 9 programs with your most common protocols
- Full-size, backlit screen
- Quick in-lab calibration for accuracy and precision
- Fully autoclavable tip cone ensures sterility and avoids cross-contamination
- Long-life battery
- Extremely lightweight construction
- Trigger action pipetting - index finger controls dispensing with just a soft press of a button
- 120° adjustable finger rest for balance and comfort
- Soft-touch tip ejection minimises thumb strain; can be operated by either hand



Code	Alt Ref	Channel	Volume, µL	Colour Code	Price
PIP1700	46200000	Single	1 - 10	Pink	£800.00
PIP1702	46200100	Single	1 - 10	Yellow	£800.00
PIP1704	46200200	Single	5 - 50	Turquoise	£800.00
PIP1706	46200300	Single	5 - 50	Yellow	£800.00
PIP1708	46200400	Single	10 - 100	Yellow	£800.00
PIP1710	46200500	Single	30 - 300	Orange	£800.00
PIP1712	46200600	Single	100 - 1000	Blue	£800.00
PIP1714	46200700	Single	500 - 5000	Green	£800.00
PIP1716	46200800	Single	100 - 10000	Green	£800.00
PIP1718	46300000	8	1 - 10	Pink	£1280.00
PIP1720	46300200	8	5 - 50	Yellow	£1280.00
PIP1722	46300400	8	30 - 300	Orange	£1258.00
PIP1724	46300800	8	100 - 1200	Turquoise	£1302.00
PIP1726	46300100	12	1 - 10	Pink	£1449.00
PIP1728	46300300	12	5 - 50	Yellow	£1449.00
PIP1730	46300500	12	30 - 300	Orange	£1449.00
PIP1732	46300700	16	5 - 50	Turquoise	£1710.00

Accessories

Code	Description	Capacity	Price
PIP1734	Carousel stand	3 Novus pipettes	£206.00

