

Gilson MICROMAN® E Pipettes Continued



Capillary Pistons (CPs)

Designed to fit the Gilson MICROMAN® positive displacement pipette, capillary piston (CP) tips are ideal for difficult to pipette samples including; viscous, volatile, dense, hot or cold, or hazardous liquids. The tips work like a syringe providing protection against pipette, sample, and user contamination, and produce accurate, reproducible, and precise results. Made of pure polypropylene, their slim fit allows you to use even the narrowest and deepest test tubes.

- Maximum accuracy and precision
- No contact between the sample and the pipette protects against contamination and provides protection for the user and the pipette, even with aerosols or corrosive vapors
- Compact and convenient TIPACK™ packaging contains colour coded racks for easy identification and compliance with Good Laboratory Practice

Code	Description	Volume Range, µL	Sterility	Pack	Price
F148412	CP10	1 - 10	Non sterile	192 (2 x TIPACKs™ of 96)	£72.00
F148312	CP10	1 - 10	Non sterile	960 (10 x TIPACKs™ of 96)	£280.00
F148413	CP10ST	1 - 10	Sterile	192 (2 x TIPACKs™ of 96)	£84.00
F148313	CP10ST	1 - 10	Sterile	960 (10 x TIPACKs™ of 96)	£408.00
F148012	CP25	3 - 25	Non sterile	960 (10 x TIPACKs™ of 96)	£306.00
F148112	CP25	3 - 25	Non sterile	200 (bulk)	£60.00
F148712	CP25ST	3 - 25	Sterile	576 (6 x TIPACKs™ of 96)	£245.00
F148113	CP50	20 - 50	Non sterile	200 (bulk)	£60.00
F148013	CP50	20 - 50	Non sterile	960 (10 x TIPACKs™ of 96)	£306.00
F148713	CP50ST	20 - 50	Sterile	576 (6 x TIPACKs™ of 96)	£245.00
F148414	CP100	10 - 100	Non sterile	192 (2 x TIPACKs™ of 96)	£72.00
F148314	CP100	10 - 100	Non sterile	960 (10 x TIPACKs™ of 96)	£280.00
F148415	CP100ST	10 - 100	Sterile	192 (2 x TIPACKs™ of 96)	£84.00
F148315	CP100ST	10 - 100	Sterile	960 (10 x TIPACKs™ of 96)	£408.00
F148114	CP250	50 - 250	Non sterile	200 (bulk)	£60.00
F148014	CP250	50 - 250	Non sterile	960 (10 x TIPACKs™ of 96)	£306.00
F148714	CP250ST	50 - 250	Sterile	576 (6 x TIPACKs™ of 96)	£245.00
F148560	CP1000	100 - 1000	Non sterile	182 (2 x TIPACKs™ of 91)	£72.00
F148180	CP1000ST	100 - 1000	Sterile	182 (2 x TIPACKs™ of 91)	£84.00



Gilson DISTRIMAN® Positive Displacement Repetitive Pipette



A positive-displacement, continuously adjustable, repetitive pipette with direct read out designed to simplify multiple dispensing when used with DISTRITIPS® syringes. With DISTRIMAN®, you can dispense the exact volume required for any protocol, even fractional volumes.

- Dispense exact volumes without tedious calculations
- Volumes are displayed clearly in microliters or milliliters on the digital read out
- Dispense the exact volume required for any analysis protocol and achieve optimal pipetting results even with volatile or viscous liquids
- No risk of accidental conversion errors with the direct read out
- Precise and reliable performance

Code	Volume Range	Price
F164001	1µL - 1.25mL	£402.00