



Corning® Step-R™ Repeating Pipettors

CORNING


Suitable to dispense volatile, dense, viscous, or foaming liquids, and for all applications where long dosing series are required. It works with 8 syringe tips sizes to dispense doses from 1µL up to 5mL.

- Easy volume setting - volume selection chart conveniently placed on the handle
- Fast and precise multiple dispensing of up to 48 pipetting doses without refilling the syringe tip
- Maintenance free
- Each device is tested for functionality and precision according to EN ISO 8655
- Supplied with an individual QC certificate

Syringe tips are available - see SLS website for further details

Code	Description	Supplied with	Price
4217	Step-R™ repeating pipettor	25 and 50mL syringe tip adapters	£356.00
6636	Step-R™ starter pack	25 and 50mL syringe tip adapters, syringe tips (5 x 0.05mL, 5 x 2.5mL, 5 x 5mL, 5 x 12.5mL, 2 x 25mL, 2 x 50mL)	£434.00

SLS Select Pipette Tips



Bulk Tips


- Packaged in resealable bags
- Suitable for use with most brands of pipette

Supplied in packs of 1000

Code	Description	Volume, µL	Price
SLS5039	White	10	£34.70
PIP7812	Yellow, basic	200	£17.40
SLS5038	Yellow, premium	200	£20.30
PIP7806	Yellow, bevelled, premium	200	£19.00
PIP7880	Blue, basic	1000	£20.80
PIP7836	Blue, bevelled	1000	£24.00
SLS5046	Blue, premium	1000	£29.90

Racked Tips

- Hinged lid
- Solid platform-great for loading tips on multichannel pipettes
- Racks can be autoclaved
- Individually housed tip further reduces contamination and tip-to-tip static attraction



A

Code	Description	Volume, µL	Sterility	Pack	Price
SLS5032	Long reach, low binding	10	Non sterile	10 x 96	£54.00
SLS5034	Long reach, low binding	10	Sterile	10 x 96	£55.00
SLS5040	Low binding	200	Non sterile	10 x 96	£51.00
SLS5042	Low binding	200	Sterile	10 x 96	£54.00
SLS5048	Long reach	1000	Non sterile	8 x 96	£52.00
SLS5050	Long reach	1000	Sterile	8 x 96	£59.00