

Corning® DeckWorks™ Pipette Tips

CORNING



A

Universal Fit Tips

Racked tips available in 10, 200, 300 and 1000µL sizes. The 10, 300 and 1000µL tips are extended length and allow for improved access into micro (1.5mL) and standard conical bottom (15/50mL) centrifuge tubes. Robust 96 tip racks are designed for repeated steam autoclave cycles and can be reused by removing the empty grey tip deck and replacing it with a new reload deck of 96 tips from the DeckWorks™ reload system. Racked 300µL DeckWorks™ tips are designed for all popular 300µL multichannel instruments.

Code	Volume Range, µL	Sterility	Pack	Price
4115	0.1 - 10	Non sterile	40 x 96	£233.00
4120	0.1 - 10	Sterile	40 x 96	£255.00
4116	1 - 200	Non sterile	40 x 96	£233.00
4121	1 - 200	Sterile	40 x 96	£255.00
4118	1 - 300	Non sterile	40 x 96	£267.00
4123	1 - 300	Sterile	40 x 96	£293.00
4119	100 - 1000	Non sterile	768 x 4	£267.00
4124	100 - 1000	Sterile	768 x 4	£293.00



A

Low Binding Filter Tips

Feature an inert, high-density polyethylene sterile filter.

Code	Volume Range, µL	Pack	Price
4135	0.1 - 10	40 x 96	£459.00
4136	1 - 20	40 x 96	£459.00
4137	1 - 100	40 x 96	£459.00
4139	1 - 200	32 x 96	£459.00
4138	1 - 200	40 x 96	£459.00
4140	100 - 1000	32 x 96	£459.00



A



Reload Systems

Available in 10, 200, 300 and 1000µL sizes and allow for convenient and economical reloading of Corning® DeckWorks™ hinged racks. Packaged in transparent trays for easy tip identification, require no transfer devices, save space and minimise packaging waste. Reload decks and rack bases are manufactured from 100% recycled polypropylene. Non sterile decks and racks are steam autoclavable.

Code	Description	Volume Range, µL	Pack	Price
4125	Reload	0.1 - 10	40 x 96	£169.00
4126	Reload	1 - 200	40 x 96	£169.00
4128	Reload	1 - 300	40 x 96	£169.00
4129	Reload	100 - 1000	32 x 96	£169.00
4153	Reload, low binding	0.1 - 10	40 x 96	£244.00
4154	Reload, low binding	1 - 200	40 x 96	£244.00
4155	Reload, low binding	1 - 300	40 x 96	£244.00