



## Rodwell Phoenix Autoclaves



The Phoenix is the ideal compact autoclave for any small scale operation that needs reliable functionality without losing valuable space.

- Front loading, cylindrical chamber
- MP25 controller
- Robust space-saving design
- Single phase "plug & play"
- Fits on a standard bench
- Integral drip tray
- Grade 316 stainless steel chamber
- Anti-door bolt over tight mechanism
- Complete versatility for almost any application

Code	Alt Ref	Capacity, L	Dims, w x d x h, mm	Chamber Dims, Ø x d, mm	Weight, kg	Media Bottle Capacity	Price
AUT1290	PHX40	40	550 x 710 x 790	350 x 460	108	12 x 1L/4 x 2L	£9625.00
AUT2060	PHX60	60	550 x 860 x 790	350 x 600	115	17 x 1L/6 x 2L	£9900.00

### Options and Accessories

Code	Alt Ref	Description	Price
AUT2000	AC27	Waste discard system	£818.00
AUT1357	AC4	Fan assisted cooling	£464.00
AUT2019	AC164	Condense collector	£277.00
AUT1374	AC276	Shelf support plus 3 additional shelves	£751.00



## THERMAL WASTEWATER STERILISATION

### The AstellBio Sink

Liquids washed down the sink are sterilised with heat before being despatched to the drain

### The AstellBio Micro EDS

A plumb-in effluent decontamination system which automatically sterilises wastewater using heat



### The AstellBio Sink & Autoclave Combo

The perfect sterilisation station for autoclaveable materials and liquid waste

AstellBio manufacture thermal effluent decontamination systems (EDS) for sterilising up to 6300L per day. Visit [AstellBio.com](http://AstellBio.com) to find out more.

[AstellBio.com](http://AstellBio.com)

**Astell**  **Bio**