

## Nickel Electro Clifton® NE5 Shaking Water Baths Continued



## Flask Clips

28L bath requires separate purchase of accessory tray to fit flask clips.

Code	Alt Ref	For use with	Price
<b>BAT7242</b>	FC0025	25mL flask	<b>£16.00</b>
<b>BAT7244</b>	FC0050	50mL flask	<b>£17.00</b>
<b>BAT7246</b>	FC0100	100mL flask	<b>£17.00</b>
<b>BAT7248</b>	FC0250	250mL flask	<b>£24.00</b>
<b>BAT7250</b>	FC0500	500mL flask	<b>£27.00</b>
<b>BAT7256</b>	FC3000	3L flask	<b>£34.90</b>
<b>BAT7259</b>	M11	Microplate	<b>£40.10</b>



## Trays

The universal spring clip tray is a multi-spring tray for accommodating a variety of flasks and bottle types, even different sizes. Featuring a criss cross pattern of springs flasks and bottles can be easily inserted and removed. 28L baths are supplied with a trolley platform to accommodate universal spring tray, flask clip trays and test tube rack trays. The tray simply drops onto the reciprocating trolley and locates snugly.

Code	Alt Ref	Description	Price
<b>BAT7177</b>	SFT9000	Universal spring tray	<b>£228.00</b>
<b>BAT7240</b>	SBF01804	Perforated accessory tray	<b>£106.00</b>

## Grant LT Ecocool Recirculating Baths



Innovative, eco-friendly, refrigerated heated circulating baths offering significant running cost savings whilst delivering powerful cooling. Consisting of two models, the LT ecocool range is supplied assembled as ready to use kits, complete with accessory hosing, clips and connectors as standard.

- Active cooling through the whole temperature range
- True energy saving of up to 80% against standard compressor units
- Thermostat and chiller work in harmony, neither will operate alone, eliminating any danger of overheating or freezing
- Single front switch for user convenience
- Modern, sleek, attractive design
- Supplied assembled as ready to use kits, complete with accessory hosing, clips and connectors as standard



Code	Alt Ref	Capacity, L	Max Flow Rate, L/min	Max Pump Pressure, mbar	Display	Temp Range, °C	Price
<b>BAT3740</b>	LT ecocool 100	5	17	250	4 digit LED	-20 to 100	<b>£2999.00</b>
<b>BAT3742</b>	LT ecocool 150	6	14 - 22 (adjustable)	530	Full colour QVGA TFT	-25 to 150	<b>£3370.00</b>