



## Sartorius Optifit Pipette Tips Continued

SARTORIUS

### Standard Tips in Refill Towers

Code	Alt Ref	Volume Range, $\mu\text{L}$	Length, mm	Pack	Price
PIP2054	790012	0.1 - 10	31.5	10 x 96	£50.00
PIP2066	790352	5 - 350	54	10 x 96	£44.00
PIP2072	791002	10 - 1000	71.5	10 x 96	£53.00
PIP2082	791202	50 - 1200	71.5	10 x 96	£54.00

### Standard Tips in Flexibulk Packs

- Combines the low cost of bulk tips with saving space, improving usability and securing tip purity
- Reduces up to 40% on pack volume
- Cuts packaging waste by nearly 50%
- Provides quick "pick & place" convenience for loading tips into tip racks

Code	Alt Ref	Volume Range, $\mu\text{L}$	Length, mm	Pack	Price
PIP2040	LH-B790204	0.5 - 200	51	960	£28.00
PIP2042	LH-B790354	5 - 350	54	960	£30.00
PIP2044	LH-B791004	10 - 1000	71.5	480	£16.00
PIP2046	LH-B791024	10 - 1000*	68.5	480	£20.00
PIP2048	LH-B791204	50 - 1200	71.5	480	£19.00

\* Wide bore tip



## Sartorius Large Volume Pipette Tips

SARTORIUS

- Certified free from RNase, DNase, DNA and PCR inhibitors
- Autoclavable
- Supplied in bags

Code	Alt Ref	Volume Range, mL	Pack	Price
PIP7736	780300	5	100	£14.00
PIP7737	LH-780316	10	250	£50.00



## Sartorius SafetySpace® Filter Tips

SARTORIUS

SafetySpace® filter tips produce repeatable results with Sartorius pipettes and eliminate the risk of your sample ever touching the filter. The unique SafetySpace® between the sample and the filter enables you to use the full volume of the tip with any pipetting mode.

- Made of virgin polypropylene (tip) and polyethylene (filter)
- Filter barrier helps protect against cross-contamination of the sample and/or the pipette and reduces the need for pipette maintenance
- Packed in color coded single tray racks protected by an airtight wrapping
- DNase, RNase, human DNA and endotoxin free certified
- Pre-sterilised

Single tray packages, wrapped in air-tight plastic

