

Corning® Reagent Reservoirs



- Manufactured from modified polystyrene
- White
- Sterile
- Disposable

Code	Capacity, mL	Inner Pack	Pack	Price
4871	50	1	100	£279.00
4870	50	5	200	£241.00
4873	100	1	100	£329.00
4872	100	5	200	£288.00
4877	12 channel (approx 5mL)	1	100	£195.00
4878	8 channel (approx 7mL)	1	100	£190.00



Blue-Ray Biotech EzMate™ Automatic Pipetting Systems



Ideal for routine low volume pipetting. It helps eliminate tedious pipetting and human error, which enables more precise and standardised results to improve productivity. With highly-customised adapters and application-oriented software, EzMate™ can satisfy the various needs of assays, including PCR/qPCR set up, sample normalisation, NGS preparation, and reagent kit manufacturing methods.

- Easy to use system setup software without need of additional technical support
- Built in automatic PCD/qPCR set-up protocols that support USB data transfer
- Single, 4 and 8 channel, 50µL or 200µL automated pipetting modules (APMs)
- 6 interchangeable standard microplate/tip rack adapters and 2 interchangeable reagent vessel adapters
- Active cooling and heating module (ACHM) options for sensitive samples/reagents
- UV lamp and HEPA filter meet the need of high-contamination assays
- Customised adapters support various types of labware
- Excellent results for qPCR standard curve and replicates
- Better precision and pipetting results than manual pipetting
- A wide variety of application-oriented method support, including PCR/qPCR set up, NGS preparation, sample normalisation, reagent kit manufacturing, cell-based assay, etc.
- Supplied with 15.6" notebook computer and EzStarter control software



Code	Alt Ref	Description	Dims, w x d x h, mm	Price
PIP1922	AREZ401-100000	EzMate 401 4-position	590 x 450 x 470	£11242.00
PIP1924	EZ4SU00-100000	EzMate 401s 4-position with UV and HEPA	660 x 450 x 490	£12795.00
PIP1926	EZ4SUA0-100000	EzMate 401s 4-position with UV and HEPA and 2-position active cooling/heating module	660 x 450 x 490	£17764.00
PIP1928	EZ4SUAC-100000	EzMate 401s 4-position with UV and HEPA and 3-position active cooling/heating module	660 x 450 x 490	£19627.00
PIP1930	AREZ601-100000	EzMate 601s 6-position	740 x 450 x 470	£13416.00
PIP1932	EZ6SU00-100000	EzMate 601s 6-position with UV and HEPA	810 x 450 x 490	£15280.00
PIP1934	EZ6SUA0-100000	EzMate 601s 6-position with UV and HEPA and 2-position active cooling/heating module	810 x 450 x 490	£20248.00
PIP1894	EZ6SUAC-100000	EzMate 601s 6-position with UV and HEPA and 3-position active cooling/heating module	810 x 450 x 490	£22112.00

Accessories

Please enquire for full range of accessories.

Code	Alt Ref	Description	Price
PIP1896	AREZ401-100000	1 channel 50µL pipetting module (new bracket)	£932.00
PIP1902	EZ4SU00-100000	1 channel 200µL pipetting module (new bracket)	£932.00
PIP1908	EZ4SUA0-100000	4 channel 50µL pipetting module (new bracket)	£931.00
PIP1906	EZ4SUAC-100000	4 channel 200µL pipetting module (new bracket)	£1180.00
PIP1898	EZ6SUAC-100000	8 channel 50µL pipetting module (new bracket)	£1180.00
PIP1900	AREZ601-100000	8 channel 50µL pipetting module for Axygen FX-50-R (new bracket)	£1180.00
PIP1904	EZ6SU00-100000	8 channel 200µL pipetting module (new bracket)	£1180.00