




## SLS Lab Pro Digital Ultrasonic Baths



With capacities from 4.5 to 25L, the SLS Lab Pro digital ultrasonic baths offer precise and accurate process control of time, temperature and ultrasonic power through the LCD control panel. Incorporating frequency LEAP technology for homogeneous, repeatable ultrasonic activity throughout the tank. The SD card port means that each cycle can be validated for full process traceability via transfer to computer and provided software. Supplied with a stainless steel basket and plastic lid.

- Ambient to 70°C temperature control
- 0 – 99 minute timer
- Degas mode
- Energy saving mode
- Consistent cleaning results every time

Code	Capacity, L	Overall Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
<b>SLS2030</b>	4.5	345 x 200 x 265	300 x 150 x 150	<b>£1062.00</b>
<b>SLS2038</b>	9.5	550 x 185 x 265	500 x 140 x 150	<b>£1318.00</b>
<b>SLS2040</b>	12.5	345 x 285 x 365	295 x 245 x 200	<b>£1515.00</b>
<b>SLS2042</b>	13	375 x 345 x 265	330 x 300 x 150	<b>£1527.00</b>
<b>SLS2044</b>	17.5	375 x 345 x 365	330 x 300 x 200	<b>£1672.00</b>
<b>SLS2046</b>	25	550 x 345 x 365	500 x 300 x 200	<b>£1965.00</b>



## Polypropylene Floating Spheres

An ideal alternative to a lid for reducing thermal losses from open tanks of water. A blanket of hollow spheres covers 91% of the liquid surface, trapping air in the gaps which adds to the thermal insulation properties and reduces steaming with no impedance to workload.

Code	Pack	Price
<b>BAT3110</b>	300 spheres	<b>£74.00</b>



## Nickel Electro Clifton® Bath Armor™



Tackles water contamination by utilising a dry, metallic thermal media comprised of small beads instead of water. Thermally efficient pellets that are easy to keep clean, Bath Armor™ eliminates nearly all water bath cleaning, protects the bath from drying out and avoids the need for racks, floaties and bottle weights. Can be cleaned using a disinfectant solution and or autoclaved.

Code	Alt Ref	Pack	Price
<b>BAT1005</b>	BA-5	4L	<b>£440.78</b>