



Gilson MINIPULS® 3

GILSON®

A high performance, low pulse, peristaltic pump with interchangeable pump heads for delivering a diverse variety of liquids - from biological to chemical - without altering the critical components within the liquid. Utilising positive displacement pump technology, liquid flows nearly pulse free as a result of as many as ten stainless steel rollers. The MINIPULS® 3 produces smoother flows with excellent accuracy and precision, even at high flow rates.

- Choice of quickly interchangeable pump heads (sold separately)
- Easy to install, simple to operate, and easy to maintain
- Uncomplicated keypad control and easy to read LCD screen
- Remote control is also possible
- Adjustable speed stepper motor and high quality precision bore tubing assure reproducible flow rates
- Supplied as drive unit only - pump heads required sold separately

Code	Price
F155001	£1487.00

Pump Heads

Code	Description	Max Flow Rate, mL/min	Channel	Price
F117604	R1 medium flow pump head	30	1	£618.00
F117800	R2 medium flow pump head	30	2	£712.00
F117606	R4 medium flow pump head	30	4	£985.00
F117608	R8 medium flow pump head	30	8	£1498.00
F117831	R4 high flow pump head	220	4	£1570.00

IT'S IN OUR DNA...



0115 982 1111



www.scientificlabs.co.uk



**SCIENTIFIC
LABORATORY
SUPPLIES**