

## Nickel Electro Clifton® Test Tube Racks



Stainless steel 304 spec 18/8 non magnetic with bright finish. Supplied flat packed ready for assembly. Can be easily disassembled for cleaning and sterilisation purposes.

Code	Alt Ref	No. of Places	Hole Ø, mm	Price
<b>BAT7187</b>	6900	12	32	£52.00
<b>BAT7186</b>	6871	16	26	£52.00
<b>BAT7184</b>	6873	18	19	£52.00
<b>BAT7182</b>	6870	26	17	£52.00
<b>BAT7180</b>	6872	36	13	£52.00



## JULABO Aqua Stabil Protective Water Bath Media



Stops the growth of algae, bacteria and other microorganisms with its bactericidal effectiveness. It protects the bath fluid and provides a long period of use for the fluid.

- Add 2mL Aqua Stabil per litre of water
- Whenever the blue tinting disappears, but latest after three months, fill the bath with fresh water and Aqua Stabil
- Batch and best before date on the bottom of the bottle

Code	Alt Ref	Pack	Price
<b>BAT5584</b>	8940006	6 x 100mL	£209.00
<b>BAT5586</b>	8940012	12 x 100mL	£356.00



## Thermo Scientific™ Aqua-Tec Water Preservation Cells



Proven to inhibit the growth of a wide range of commonly encountered microbes for up to six months. Eliminate the need to handle dangerous chemicals or expose sensitive cells to the risk of potentially toxic off-gassing from disinfectants. The Aqua-Tec cell allows your CO<sub>2</sub> incubators and water baths to operate at maximum efficiency and saves the significant time, handling, and exposure to chemicals that come with using traditional chemical disinfection methods.

- Designed and tested specifically for laboratory use in CO<sub>2</sub> incubators and water baths
- Proven effective for up to 6 months against bacteria, moulds, and fungi
- Safe to handle – requires no special disposal protocol
- Self regulates to maintain effectiveness, regardless of water level
- Easy to use in equipment from any manufacturer

Code	Pack	Price
<b>WAT1200</b>	4	£363.00

