



Tubing Adapters, Right Angle



- Allows the connection of flexible laboratory tubing to a jointed glassware assembly
- With right angle connection for attachment of rubber hose etc.
- Detachable plastic screwthread connectors allow flexible tubing to be fitted easily and safely

Code	Cone Size	Price
QMF102B	19/26	£19.71
QMF103B	24/29	£19.86
QMF104B	29/32	£24.40
QMF105B	34/35	£28.60



Cone/Screwthread Adapters



- With screwthread joint complete with plastic cap, silicone rubber ring and PTFE washer
- Ideal for securing thermometers, gas inlet tubes, etc.

Supplied in packs of 2

Code	Cone Size	Thread	Ø Through Cap, mm	Price
QST51130	14/23	13	6.0 - 7.0	£40.00
QST51180	14/23	18	7.0 - 8.5	£40.00
QST52130	19/26	13	6.0 - 7.0	£40.00
QST52180	19/26	18	7.0 - 8.5	£40.00
QST53130	24/29	13	6.0 - 7.0	£49.50
QST53180	24/29	18	7.0 - 8.5	£49.50




Stopcock Adapters with Cone Joint



- Straight pattern adapters
- Ground conical cone joint on one end and a detachable plastic screwthread connector on the other
- Rotaflo® General Purpose PTFE stopcock to control the flow of liquids and gases (bore: 3mm)

Code	Cone Size	Price
QMF111	14/23	£61.00
QMF112	19/26	£57.00
QMF113	24/29	£57.00



T Shape Stopcock Adapters with Cone Joint



- T shape stopcock adapters
- Ground conical cone joint on one end and a detachable plastic screwthread connector on the other
- Rotaflo® General Purpose PTFE stopcock to control the flow of liquids and gases (bore: 3mm)

Code	Socket Size	Price
QMF111BSC	14/23	£69.00
QMF112BSC	19/26	£60.00
QMF113BSC	24/29	£61.00