

Straight Receiver Adapters



- Used to safely feed the distillate from condensers into the receiving vessels
- Delivery receiver with straight stem
- Approx overall height: 180mm

Code	Socket Size	Price
QRA122	19/26	£21.46
QRA123	24/29	£21.46



Plain Bend Socket to Cone Receiver Adapters



- Without vent
- Bend angle: 100°

Code	Socket Size	Cone Size	Price
QRA911	14/23	14/23	£29.80
QRA922	19/26	19/26	£29.60
QRA933	24/29	24/29	£29.60
QRA944	29/32	29/32	£35.50



Water Estimators, Light Entrainer Receivers



- Complies with BS 756
- Dean & Stark pattern receivers for the determination of the water content in petroleum products
- For use with solvent entrainers which are lighter than water
- Borosilicate glass with graduated tube
- Receiver QWE62310T is fitted with PTFE GP Rotaflo® stopcock for easy drainage
- Cone size: 24/29, socket size: 19/26

Code	Description	Capacity, mL	Graduations, mL	Price
QWE6232	Receiver	2.0	0.05	£73.00
QWE6237	Receiver	7.5	0.10	£82.00
QWE62325	Receiver	25.0	0.20	£102.00
QWE62310T	Receiver with stopcock	10.0	0.10	£144.00



Water Estimators, Heavy Entrainer Receivers



- Complies with BS 756
- For the determination of water content in petroleum products using a solvent entrainer which is heavier than water
- Also suitable for the estimation of volatile oil content
- Highest quality borosilicate glass with graduated tube and QUICKFIT® precision-ground joints
- Cone size: 24/29, socket size: 19/26

Code	Capacity, mL	Graduations, mL	Price
QWE7233	3.0	0.05	£113.00
QWE72312	12.5	0.10	£110.00

