

## Soxhlet Extractors



- Conforms to BS 2071
- For extracting soluble components or impurities from solids via repeated soaking/siphoning with solvent
- Accurately formed concentric siphon tube for reliable performance

Code	Nominal Capacity (Siphoning Volume), mL	Socket Size	Cone Size	Recommended Thimble Size, mm	Price
<b>QEX555100</b>	100	34/35	34/35	28 x 120	<b>£91.00</b>
<b>QEX563</b>	100	40/38	24/29	33 x 100	<b>£88.00</b>
<b>QEX583</b>	200	50/42	24/29	41 x 123	<b>£109.00</b>



## Round Bottom Short Neck Flasks



- With highest quality QUICKFIT® ground glass socket
- Manufactured from borosilicate glass for corrosion-free performance and resistance to thermal shock
- Body diameter and wall thickness comply with relevant ISO 1773 and ISO 4797 requirements

Code	Nominal Capacity, mL	Socket Size	Approx Overall Height, mm	Approx Neck Length, mm	Pack	Price
<b>QFR100150</b>	100	14/23	105	15	5	<b>£84.00</b>
<b>QFR100250</b>	100	19/26	105	15	5	<b>£84.00</b>
<b>QFR100350</b>	100	24/29	105	15	5	<b>£84.00</b>
<b>QFR250150</b>	250	14/23	126	15	2	<b>£36.80</b>
<b>QFR250250</b>	250	19/26	140	15	5	<b>£92.00</b>
<b>QFR250350</b>	250	24/29	140	15	5	<b>£92.00</b>



## Flat Bottom Short Neck Flasks



- Complies with ISO 1773 and ISO 4797
- With highest quality QUICKFIT® ground glass socket
- Borosilicate glass for highest resistance to thermal shock and chemical corrosion
- Flat base allows use with hotplates
- Not suitable for use under vacuum
- Socket size: 24/29

Code	Nominal Capacity, mL	Approx Overall Height, mm	Approx Neck Length, mm	Pack	Price
<b>QFF100350</b>	100	103	15	5	<b>£88.00</b>
<b>QFF150350</b>	150	110	15	5	<b>£100.00</b>
<b>QFF250350</b>	250	130	10	5	<b>£91.00</b>
<b>QFF500350</b>	500	140	15	5	<b>£109.00</b>
<b>QFF1L3U50</b>	1000	170	25	2	<b>£60.00</b>
<b>QFF2L350</b>	2000	205	25	5	<b>£311.00</b>



## Round Bottom 2 Neck Flasks



- Vertical centre neck and an angled side neck
- Fitted with two QUICKFIT® ground glass sockets
- Side neck useful for introducing gas inlet tubes and thermometers etc. into the flask
- Borosilicate glass for highest resistance to thermal shock and chemical attack

Code	Nominal Capacity, mL	Centre Socket Size	Side Socket Size	Approx Overall Height, mm	Price
<b>QFR25151A</b>	25	14/23	14/23	77	<b>£36.30</b>
<b>QFR100352A</b>	100	24/29	19/26	110	<b>£36.30</b>
<b>QFR250352A</b>	250	24/29	19/26	130	<b>£34.80</b>

