



Thermometers with Integral Cone, Spirit Filled



- Complies with ISO 386
- Integral ground glass cone on the stem which allows direct fitting into QUICKFIT® jointed apparatus
- Filled with high visibility blue spirit

Code	Cone Size	Temp Range, °C	Price
MF312S	14/23	-10 to 250	£124.00



Plain Air/Steam Inlet Tubes



- Plain borosilicate glass tube with jet
- Designed to fit the range of QUICKFIT® thermometer/inlet adapters with size 13 thread

Code	Approx OD, mm	Price
QMF5	7	£31.50



Ground Socket Test Tubes



- High quality test tubes with a precision-ground QUICKFIT® socket
- Accepts QUICKFIT® glass stoppers
- Borosilicate glass for excellent resistance to chemical attack and thermal shock
- Stoppers sold separately

Supplied in packs of 5

Code	Socket Size	Dims, OD x h, mm	Approx Capacity, mL	Price
QMF24160	14/23	19 x 150	19	£60.00
QMF24260	19/26	23 x 150	31	£60.00
QMF2430	24/29	29 x 150	54	£60.00



Hexagonal Head Glass Stoppers



- Ideally suited for use with QUICKFIT® jointed glassware and other glassware with standard ground glass sockets
- All QUICKFIT® stoppers conform to the standard tapers as outlined in ISO 383
- Borosilicate glass for chemical durability and high temperature performance
- Hollow-blown construction to reduce weight
- Hexagonal head prevents the stopper from rolling and provides a good grip
- Top of stopper is flat, allowing them to be stood on end and avoid contamination of the ground surface

Supplied in packs of 5

Code	Cone Size	Overall Height, mm	Price
QSB100	10/19	46	£28.20
QSB120	12/21	52	£32.90
QSB140	14/23	56	£27.20
QSB190	19/26	61	£27.90
QSB240	24/29	67	£31.30
QSB290	29/32	75	£36.10