




Heathrow Scientific Coated Wire Racks



- Chemically resistant HDPE coated wire racks
- Suitable for use down to -40°C
- Autoclavable
- Two levels of grids help keep tubes/bottles secure

A

Code	Alt Ref	Bottom	Well Ø, mm	Array	No. of Places	Colour	Dims, w x l x h, mm	Price
RAC0012	HS120089	Slotted	10 - 13	6 x 12	72	Blue	105 x 205 x 65	£36.40
RAC0050	HS120759	Slotted	10 - 13	9 x 12	108	Blue	161 x 204 x 65	£31.10
RAC0052	HS120760	Slotted	13 - 16	4 x 10	40	Blue	83 x 200 x 59	£26.80
RAC0054	HS120761	Slotted	13 - 16	4 x 12	48	Blue	127 x 237 x 65	£28.90
HS23072	HS23072	Slotted	13 - 16	6 x 12	72	Blue	127 x 237 x 65	£39.70
HS23108	HS23108	Slotted	13 - 16	9 x 12	108	Red	188 x 241 x 65	£31.10
RAC0056	HS120763	Slotted	13 - 16	9 x 12	108	Blue	188 x 241 x 65	£31.10
HS23215	HS23215	Slotted	16 - 20	3 x 5	15	Blue	76 x 127 x 76	£29.80
RAC0058	HS120764	Slotted	16 - 20	4 x 10	40	Blue	90 x 240 x 80	£33.20
RAC0060	HS120765	Slotted	16 - 20	8 x 10	80	Blue	190 x 240 x 80	£34.30
RAC0068	HS120769	Slotted	20 - 25	3 x 8	24	Blue	105 x 273 x 80	£29.10
RAC0062	HS120766	Slotted	20 - 24	4 x 10	40	Blue	107 x 260 x 80	£33.20
RAC0064	HS120767	Slotted	22 - 25	4 x 10	40	Blue	118 x 290 x 80	£34.30
RAC0066	HS120768	Slotted	22 - 25	8 x 10	80	Blue	230 x 290 x 80	£35.40
HS232461	HS232461	Slotted	30	2 x 4	8	Blue	80 x 147 x 89	£29.80
RAC0070	HS120770	Slotted	30	6 x 8	48	Blue	209 x 277 x 98	£39.60
RAC0010	HS23224	Slotted	58 - 60	2 x 3	6	Blue	135 x 205 x 104	£36.40
RAC0002	HS120085	Mesh	25	2 x 6	12	Blue	64 x 164 x 32	£24.70
RAC0004	HS120086	Mesh	25	5 x 10	50	Blue	145 x 271 x 32	£41.80
RAC0006	HS120087	Mesh	33	2 x 6	12	Blue	80 x 212 x 60	£28.60
RAC0008	HS120088	Mesh	33	5 x 10	50	Blue	189 x 355 x 60	£45.00
RAC0072	HS120771	Mesh	40	3 x 8	24	Blue	139 x 357 x 98	£33.20
RAC0074	HS120772	Mesh	40	6 x 8	48	Blue	260 x 350 x 98	£37.50



Heathrow Scientific Adapt-a-Rack®



- Flexible multi tube rack where each slot adapts to accommodate 12 - 30mm (5 - 50mL) tube sizes and holds them in a firm upright position
- Racks can be linked together for additional holding
- Available in solid or a combination of colours for easy coding identification
- Made from POM (polyoxymethylene)
- Can hold different sized tubes in each slot
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Linked racks should be supported underneath when carrying
- Autoclavable
- Dims, L x w x h: 181 x 56 x 76mm

Supplied in packs of 2

A