



## Heathrow Scientific 96-Well Reversible Racks



- Easy-to-grasp tabs for sample handling and storage
- Polypropylene
- Wells are arranged in standard 8 x 12 array
- Tabs on each end of rack snap lid into place
- 0.5mL wells separated by moulded-in grid lines, making them easier to see
- Raised rings around the 1.5 to 2.0mL wells make tubes easier to remove
- Autoclavable
- Dims, L x w x h: 246 x 121 x 50mm

Supplied in packs of 5

	Side 1	Side 2
Wells	96	96
Well Ø, mm	11	7.6
For Tubes, mL	1.5/2.0	0.5

Code	Alt Ref	Colour	Price
RAC0660	HS2345B	Natural	£91.00
RAC0662	HS2345C	Blue	£91.00
RAC0664	HS2345D	Green	£91.00
RAC0666	HS2345E	Pink	£91.00
RAC0668	HS2345F	Yellow	£91.00
RAC0670	HS2345A	Assorted	£91.00



## Heathrow Scientific 80-Well Microtube Racks



- Rack to easily fill and transfer 80 tubes in compact arrangement
- Popular fraction-collector-plate style
- Holds 1.5 - 2.0mL microtubes in a 5 x 16 array
- Polypropylene
- Tubes sit in large wells with conical bottoms
- Moulded-in coordinates to aid organisation
- Autoclavable
- Dims, L x w x h: 225 x 67 x 28mm

Supplied in packs of 5

Code	Alt Ref	Colour	Price
RAC0900	HS29025A	Natural	£36.40
RAC0902	HS29025B	Blue	£36.40
RAC0904	HS29025C	Green	£36.40
RAC0924	HS29025D	Pink	£36.40
RAC0926	HS29025E	Yellow	£36.40
RAC0922	HS29025F	Orange	£36.40
RAC0917	HS29025J	Red	£36.40
RAC0916	HS29025G	Assorted	£36.40