




## Heathrow Scientific Low Temperature PCR Storage Racks



- Rack slides out from the lid, and then rotates 180°, to become a convenient stand the rack sits on
- Ergonomic textured finger grips provide a secure hold
- Lid has a side opening that enables faster cooling of samples to ensure integrity and is made from translucent polycarbonate
- A location ID of numbers and letters allows for organised sample placement
- Compact, stackable footprint saves space
- Holds 96 x 0.2mL tubes (individual or strip) in an 8 x 12 array
- Can be used for freezer storage down to -80°C
- Detachable, ambidextrous lid can be switched to accommodate left or right hand orientation
- Not recommended for autoclaving
- Dims, w x l x h: 92 x 125 x 37mm

Supplied in packs of 5

Code	Alt Ref	Colour	Price
PCR1610	HS120538	Assorted	£60.00
PCR1612	HS120539	Blue	£60.00
PCR1614	HS120540	Green	£60.00
PCR1616	HS120541	Purple	£60.00



## Azlon® Polypropylene Test Tube Racks




- Economical and sturdy test tube racks supplied in self assembly flat packs
- Autoclavable at 121°C
- Will not float in water baths
- Alpha numeric place coding
- White

Also available in other colours - see the SLS website for the full range

A

Code	Tube Ø, mm	No. of Places	Dims, w x l x h, mm	Price
STA3200	13	84	126 x 265 x 75	£13.89
STA3202	16	55	126 x 265 x 75	£15.10
STA3204	20	40	126 x 265 x 75	£13.89
STA3206	25	32	126 x 265 x 88	£15.53
STA3208	30	21	126 x 265 x 88	£15.53



## SAMCO Wooden Test Tube Racks



- Hardwood, with plastic ends
- Suitable for 6 test tubes
- Packs flat for transport and storage

Code	Alt Ref	Tube Ø, mm	Price
STA3000	L330/01	22	£15.07
STA3008	L330/21	32	£7.89