

Tube Racks

- Range of anodised aluminium racks
- Suitable for a wide range of flat or round bottom tubes

Other sizes available on request - please enquire for details



Code	Description	Tube Ø, mm	Array	Dims, w x l x h, mm	Price
STA2100	Single platform	7	5 x 10	62 x 127 x 28	£14.70
STA2102	Single platform	12	4 x 12	77 x 216 x 28	£14.70
STA2114	Single platform	17.5	5 x 10	120 x 232 x 33	£30.70
STA2116	Single platform	17.5	2 x 12	53 x 277 x 33	£15.30
STA2120	Single platform	18	5 x 10	129 x 252 x 29	£18.70
STA2124	Single platform	19.5	5 x 10	133 x 265 x 29	£28.70
STA2126	Single platform	23.5	6 x 6	173 x 177 x 29	£34.80
STA2132	Single platform	30	6 x 6	205 x 208 x 40	£36.50
STA2136	Single platform	30	2 x 6	80 x 223 x 42	£33.70
STA2202	Double platform	11	5 x 10	108 x 214 x 46	£26.70
STA2210	Double platform	14	5 x 10	107 x 210 x 60	£36.80
STA2220	Double platform	14	4 x 12	87 x 250 x 60	£43.20
STA2224	Double platform	14	12 x 12	246 x 250 x 60	£102.00
STA2226	Double platform	18	2 x 10	53 x 233 x 60	£36.60
STA2238	Double platform	18	5 x 12	120 x 277 x 60	£28.70
STA2240	Double platform	18	10 x 12	233 x 277 x 60	£76.00
STA2246	Double platform	20	4 x 10	112 x 283 x 60	£54.00
STA2248	Double platform	20	5 x 10	140 x 283 x 60	£56.00

Petri Dish Racks

- Suitable for 90mm Petri dishes
- Durable polypropylene
- Rounded corners minimising risk to users
- Open design for easy access
- Suitable for repeat autoclaving
- Clips supplied for attaching an ID card to the end pillars
- Available with bases coloured in blue, red, yellow, or green - please specify when ordering

Code	Capacity	Dims, w x l x h, mm	Price
RAC0200	30 Petri dishes	100 x 294 x 172	£43.10
RAC0202	45 Petri dishes	100 x 294 x 235	£59.00

