



BioCOLD Spark Free Refrigerator



Built to a very high specification, to meet the requirement for general laboratory storage. This upright refrigerator features 5 glass safety shelves as standard.

- Reversible doors
- Replaceable door gasket
- CFC and HCFC free refrigeration system
- Manual defrost
- Spark free



Code	Alt Ref	Capacity, L	Temp Range, °C	Dims, w x d x h, mm	Price
REF2222	BIO360FRSS	360	2 to 8	600 x 655 x 1750	£886.00



Lec Basic Laboratory Refrigerator



A basic spark free laboratory underbench refrigerator with lock, offering safe and reliable storage of volatile liquids, reagents, chemicals, poisons and samples. Includes 3 glass shelves and 1 plastic drawer.



Code	Alt Ref	Model	Capacity, L	Temp Range, °C	Dims, w x d x h, mm	Price
REF1990	444441207	LR111	130	2 to 10	545 x 554 x 845	£396.00



Lec Laboratory Refrigerators



A range of safe, reliable and energy efficient fridges for laboratories. Fridge temperature set point of +6°C and operates from 2 to 10°C.

- Spark free
- External controller
- Supplied with lock and two keys
- Audible and visual high and low temperature alarms
- Reversible doors
- Auto defrost



Code	Model	Capacity, L	Dims, w x d x h, mm	Price
REF4076	LSFSR82UK	82	502 x 540 x 660	£746.00
REF4078	LSFSR107UK	107	495 x 620 x 830	£848.00
REF4080	LSFSR158UK with access port hole	158	595 x 642 x 850	£732.00
REF4048	LR1607C upright	475	695 x 910 x 1810	£3999.00