

## Lec Pharmacy Refrigerators Continued



## Control Plus Range

- Measures your true temperature with twin probes
- Air and core temperature display
- Now with the new controller making it user friendly
- New rubber sealed port hole for easier external temperature monitoring
- Available in solid and glass doors

Code	Model	Capacity, L	Dims, w x d x h, mm	Door Style	Price
REF3902	PPGR158UK	158	595 x 642 x 850	Glass	£1641.00
REF4018	PPSR158UK	158	595 x 642 x 850	Solid	£1380.00
REF3922	PPGR273UK	273	595 x 625 x 1565	Glass	£2213.00
REF4094	PPSR273UK	273	595 x 625 x 1565	Solid	£1796.00
REF4092	PPGR353UK	353	595 x 640 x 1860	Glass	£2618.00
REF4096	PPSR353UK	353	595 x 640 x 1860	Solid	£2116.00
REF4032	PG1607C* upright	444	685 x 910 x 1688	Solid	£5221.00

\* Single temperature probe



## Liebherr Underbench Refrigerator

LIEBHERR

An energy efficient unit that easily fits on or under a counter. Designed for better cooling performance and stable temperatures. Features electronic controls for precise temperature control.

- Highly efficient dynamic cooling system ensures even temperatures and rapid cooling
- Fans stop when opening the door to save valuable energy
- Large surface evaporator for a better cold air distribution
- Digital LED temperature control panel with audible and visual high temperature and open door alarms
- Reversible, self closing door can be opened from either direction
- Features a key lock and adjustable feet in front
- Crush proof steel body with white polystyrene interior

Code	Alt Ref	Capacity, L	Temp Range, °C	Dims, w x d x h, mm	Price
REF1728	FKUv 1610	141	1 to 15	600 x 615 x 830	£859.00

## KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty