

Haier Biomedical Blood Bank Refrigerators Continued

Haier Biomedical
Intelligent Protection of Life Science

Code	Alt Ref	Volume, L	Display	Capacity	Dims, w x d x h, mm	Shelf Positions	Price
REF2310	HXC-149	149	LED	60 bags	625 x 820 x 1130	2	£7496.00
REF2312	HXC-149T	149	Touchscreen	60 bags	625 x 820 x 1130	2	£9909.00
REF2314	HXC-429	429	LED	195 bags	625 x 940 x 1830	5	£10856.00
REF2316	HXC-429T	429	Touchscreen	195 bags	625 x 940 x 1830	5	£13618.00
REF2318	HXC-629	629	LED	312 bags	765 x 940 x 1970	6	£13233.00
REF2320	HXC-629T	629	Touchscreen	312 bags	765 x 940 x 1970	6	£15657.00
REF2322	HXC-1369	1369	LED	624 bags	1545 x 940 x 1970	6	£20205.00
REF2324	HXC-1369T	1369	Touchscreen	624 bags	1545 x 940 x 1970	6	£23443.00



Haier Biomedical Pharmacy Refrigerators

Haier Biomedical
Intelligent Protection of Life Science

Small under bench pharmacy refrigerators, suitable for installation in pharmacies, pharmaceutical companies, hospitals, vaccination centres and clinics.

- Audible and visual alarms for high/low temperature, door ajar and sensor error
- Forced air cooling with optimised air distribution system designed to achieve maximum temperature uniformity and stability
- Adjustable temperature range 2 to 8°C
- Microprocessor control, LED display, solid door, and temperature adjustment with an increment of 0.1°C
- Optimised space utilisation design with multi level and adjustable shelf height to accommodate different storage requirements

Code	Alt Ref	Capacity, L	Dims, w x d x h, mm	Price
SLS1130	HYC-68	68	495 x 580 x 660	£792.00
SLS1134	HYC-118	118	597 x 635 x 810	£1333.00



Gram BIOBASIC Refrigerators

The BIOBASIC range provides cost-effective biostorage solutions for a host of applications. With forced air distribution, temperature alarms, internal and external ATEX and a user friendly interface, the BIOBASIC range offers an exceptional selection of features in a very small package.

- Temperature range: 2 to 15°C
- Acoustic door alarm, high and low temperature alarms, and a voltage free contact (NO/ NC) for remote alarms.
- Full ATEX compliance - both internal and external
- Forced air distribution system makes sure the temperature inside the cabinet remains stable and uniform
- Automatic smart defrost with re-evaporation of defrost water

Code	Model	Capacity, L	Dims, w x d x h, mm	No. of Shelves	Price
REF1160	RR210	125	595 x 640 x 831	2	£1474.00
REF1162	RR310	218	595 x 640 x 1220	3	£1200.00
REF1164	RR410	346	595 x 640 x 1876	5	£1390.00