

Labcold Spark Free Benchtop Laboratory Freezer



The Labcold spark free range has a specially manufactured spark free interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it an important part of laboratory safety. Not suitable for storing open containers of volatile substances.

- Manual defrost
- Door lock with 2 keys
- 3 shelves
- Self closing, solid door
- Temperature display
- High/low audio and visual temperature alarm, door open alarm



Code	Capacity, L	Temp Range, °C	No. of Shelves	Dims, w x d x h, mm	Price
RLVF0217	66	-18 to -25	3	450 x 510 x 740	£765.00

Labcold Spark Free Upright Laboratory Freezer



The Labcold spark free range has a specially manufactured spark free interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it an important part of laboratory safety. Not suitable for storing open containers of volatile substances.

- Manual defrost
- Door lock with 2 keys
- Self closing, solid door
- Temperature display
- High/low audio and visual temperature alarm, door open alarm
- Fan circulated
- 4 rollers



Code	Capacity, L	Temp Range, °C	No. of Shelves	Dims, w x d x h, mm	Price
RLVF1517	406	-18 to -25	6	600 x 700 x 1875	£1880.00



Quality Without Compromise

Quality, value
laboratory
equipment for
everyday use

