


## PYREX® Low Form Beakers, Heavy Duty

PYREX®

- Mechanically stronger than standard beakers
- Designed for longer life and increased user safety
- Reinforced rim and spout, less prone to breakage during washing etc
- Thicker walls and heavier bases
- Manufactured from PYREX® borosilicate glass
- White graduations and marking spot



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BEA12780	1003-150	150	57 x 89	12	£86.00
BEA12800	1003-250	250	68 x 90	12	£107.00
BEA12820	1003-400	400	77 x 110	12	£121.00
BEA12840	1003-600	600	90 x 124	6	£70.00
BEA12860	1003-1L	1000	108 x 156	6	£103.00
BEA1288	1003-2L	2000	131 x 179	Each	£31.50
BEA1289	1003-3L	3000	170 x 230	Each	£154.00
BEA1290	1003-4L	4000	160 x 250	Each	£154.00
BEA1291	1003-5L	5000	180 x 300	Each	£266.00

## Simax Low Form Beakers

SIMAX®

- Manufactured from borosilicate glass
- ISO 3819, DIN 12 331, CSN 70 4031
- Consistent and even wall thickness for longevity
- With graduations and spout



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BEA10000	155/10	10	26 x 35	10	£27.40
BEA10020	155/25	25	34 x 50	10	£22.20
BEA10040	155/50	50	42 x 60	10	£20.90
BEA10060	155/100	100	50 x 70	10	£17.50
BEA10070	155/150	150	60 x 80	10	£22.20
BEA10080	155/250	250	70 x 95	10	£19.90
BEA10100	155/400	400	80 x 110	10	£26.40
BEA10120	155/600	600	90 x 125	10	£30.10
BEA10140	155/800	800	100 x 135	10	£53.00
BEA10160	155/1000	1000	105 x 145	10	£51.10
BEA1018	155/2000	2000	130 x 185	Each	£9.00
BEA1020	155/3000	3000	150 x 210	Each	£30.60
BEA1022	155/5000	5000	170 x 270	Each	£44.40