

## Glove Plus® Ultra Nitrile Gloves



- High performance protection against a wide range of chemicals and blood borne pathogens
- Durable and strong
- Textured fingertips for enhanced grip
- Below 2mg/g glove powder content
- Enhanced processing, silicon free, ultra performance
- Suitable for contact with chemicals and blood borne pathogens
- Medium/high risk category
- Duration of use: 120+ minutes
- Length: 240mm
- PPE class III
- EN 455 parts 1 & 2
- AQL 0.65
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5
- ASTM D6978-05 permeation of chemotherapy drugs passed



Supplied in packs of 100

Code	Alt Ref	Size	Price
SAF6322	UNP32	Small	£41.00
SAF6324	UNP33	Medium	£41.00
SAF6326	UNP34	Large	£41.00
SAF6328	UNP35	X large	£41.00

## Glove Plus® Lite Nitrile Gloves



- The perfect medium weight nitrile gloves
- Superb chemical barrier properties and precision grip, they are the ideal everyday glove.
- Below 2mg/g glove powder content
- Textured fingertips and precision grip
- Medium risk category
- Suitable for contact with chemicals and blood borne pathogens
- Duration of use: 30 – 120 minutes
- Length: 240mm
- Medical grade
- EN 455 parts 1 & 2
- EU PPE class III
- AQL 0.65
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5
- ASTM D6978-05 permeation of chemotherapy drugs passed



Supplied in packs of 100

Code	Alt Ref	Size	Price
SAF5682	LN302	Small	£38.90
SAF5684	LN303	Medium	£38.90
SAF5686	LN304	Large	£38.90
SAF5688	LN305	X large	£38.90

## Glove Plus® Prime Nitrile Gloves



- The ideal glove for use over short periods of time where splash contact from aggressive chemicals is unlikely
- Barrier protection against blood borne pathogens
- Below 2mg/g glove powder content
- Increased sensitivity, fingertip textured for precision grip
- Duration of use: 0 – 30 minutes
- Suitable for contact with chemicals and blood borne pathogens
- Low risk category
- Medical grade
- Length: 240mm
- EN 455 parts 1 & 2
- EU PPE class III
- AQL 0.65
- ISO 16604b:2003 viral test passed
- EN 420 dexterity test passed level 5



Supplied in packs of 100

Code	Alt Ref	Size	Price
SAF5520	PN301	X small	£18.90
SAF5522	PN302	Small	£18.90
SAF5524	PN303	Medium	£18.90
SAF5526	PN304	Large	£18.90
SAF5529	PN305	X large	£18.90