

SHIELDskin™ ORANGE NITRILE™ 260 Gloves



- Orange nitrile/neoprene
- Powder, latex and accelerator free
- Textured fingertips for grip
- Beaded cuff
- Ambidextrous
- twinSHIELD™ double walled protection
- Length: 260mm
- Category III
- AQL 0.25
- EN ISO 374-1:2016 + A1:2018 Type B JKPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 455-1/2/3/4:2000 (medical gloves for single use)



Supplied in packs of 90

Code	Alt Ref	Size	Price
SAF4310	676231	X small (6)	£59.00
SAF4312	676232	Small (7)	£59.00
SAF4314	676233	Medium (8)	£59.00
SAF4316	676234	Large (9)	£59.00
SAF4318	676235	X large (10)	£59.00

SHIELDskin™ ORANGE NITRILE™ 300 Gloves



- Orange nitrile/neoprene
- Powder, latex and accelerator free
- Textured fingertips for grip
- Beaded cuff
- Ambidextrous
- twinSHIELD™ double walled protection
- Length: 300mm
- Category III
- AQL 0.25
- EN ISO 374-1:2016 + A1:2018 Type B JKPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 455-1/2/3/4:2000 (medical gloves for single use)

Sterile

- Sterile (Sterility Assurance Level 10^{-6}) to ISO 11137-2:2015

Supplied in cases of 160 pairs (1 pair/pack, 20 packs/bag, 8 bags/case)

Code	Alt Ref	Size	Price
SAF10800	676351	X small (6)	£531.10
SAF10820	676352	Small (7)	£531.10
SAF10840	676353	Medium (8)	£531.10
SAF10860	676354	Large (9)	£531.10
SAF10880	676355	X large (10)	£473.00



Non Sterile

Supplied in packs of 50 pairs

Code	Alt Ref	Size	Price
SAF4300	676251	X small (6)	£50.00
SAF4302	676252	Small (7)	£50.00
SAF4304	676253	Medium (8)	£50.00
SAF4306	676254	Large (9)	£50.00
SAF4308	676255	X large (10)	£50.00

