

## Kimtech™ Purple Nitrile™ Xtra™ Gloves

KIMTECH™

- Purple synthetic nitrile polymer
- Powder, latex and silicone free
- Textured fingertips for enhanced grip
- Beaded cuff, ambidextrous, anti static
- Length: 300mm
- Category III
- AQL 0.65
- EN ISO 374-1:2016 Type B JKT (chemical risks)
- EN 374-4:2014 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- Food contact approved

Supplied in packs of 50



Code	Alt Ref	Size	Price
SAF3338	97610	X small	£45.00
SAF3340	97611	Small	£45.00
SAF3342	97612	Medium	£45.00
SAF3344	97613	Large	£45.00
SAF3346	97614	X large	£45.00

## Kimtech™ Prizm™ Neoprene-Nitrile Gloves

KIMTECH™

- Multi-layered, dark violet/magenta neoprene-nitrile, powder and latex free, accelerator free donning layer
- Ultra fingertip grip, beaded cuff, ambidextrous, certified ergonomic comfort and best in class dexterity
- Length: 240mm (XS - M), 250mm (L - XL)
- Category III
- AQL 0.65
- EN ISO 374-1:2016 Type A JKLMPT (chemical risks)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ASTM D6978-05 (cytotoxic drugs)



Code	Alt Ref	Size	Pack	Price
SAF1672	99221	X small	100	£54.00
SAF1674	99222	Small	100	£54.00
SAF1676	99223	Medium	100	£54.00
SAF1678	99224	Large	100	£54.00
SAF1680	99225	X large	90	£54.00

## Kimtech™ Prizm™ Xtra™ Neoprene-Nitrile Gloves

KIMTECH™

- Multi-layered, dark violet/magenta neoprene-nitrile, powder and latex free, accelerator free donning layer
- Ultra fingertip grip, beaded cuff, ambidextrous, certified ergonomic comfort and best in class dexterity
- Length: 300mm
- Category III
- AQL 0.65
- EN ISO 374-1:2016 Type A JKLMPT (chemical risks)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ASTM D6978-05 (cytotoxic drugs)

Supplied in packs of 50



Code	Alt Ref	Size	Price
SAF1682	99251	X small	£41.00
SAF1684	99252	Small	£41.00
SAF1686	99253	Medium	£41.00
SAF1688	99254	Large	£41.00
SAF1690	99255	X large	£41.00