

DURAN® Silicone Lids Continued



Sets


Sets of 1 x small, 1 x medium, and 1 x large lids.

Code	Alt Ref	Colour	Price
BOT1242	291100011	Pink	£5.80
BOT1244	291100028	Cyan	£5.80
BOT1246	291100036	Green	£5.80

Academy Low Form Beakers

ACADEMY

- Made from borosilicate glass 3.3
- Superior chemical and heat resistance
- White graduations and marking area




Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BEA34020	A/2218/50	50	42 x 60	12	£12.50
BEA34040	A/2218/100	100	50 x 70	12	£11.00
BEA34060	A/2218/250	250	70 x 95	12	£12.50
BEA34080	A/2218/400	400	80 x 110	12	£19.10
BEA34100	A/2218/600	600	90 x 125	6	£12.00
BEA3412	A/2218/800	800	100 x 135	6	£17.20
BEA3414	A/2218/1000	1000	105 x 145	6	£17.80
BEA3416	A/2218/2000	2000	130 x 190	Each	£9.90

Academy Low Form Beakers, Heavy Duty

ACADEMY

- Robust but economical heavy duty beakers
- Borosilicate 3.3 glass
- Increased wall thickness



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BEA0120	A/2283/150	150	64 x 82	12	£35.10
BEA0122	A/2283/250	250	74 x 94	12	£41.80
BEA0124	A/2283/400	400	84 x 112	12	£65.00
BEA0126	A/2283/600	600	94 x 127	6	£41.10
BEA0130	A/2283/1000	1000	109 x 147	6	£50.00
BEA0132	A/2283/2000	2000	134 x 187	4	£47.80