



## Glove Plus® Polymer Coated Latex Gloves



- High quality and cost effective
- Barrier against blood-borne pathogens
- Internal polymer coating enables easy donning and extra comfort
- Fully textured finish for additional grip
- Protein content 50µg/g glove
- Below 2mg/g glove powder content
- Ambidextrous and powder free
- Duration of use: 120+ minutes
- Suitable for use with blood borne pathogens (not suitable for contact with chemicals)
- Medium/high risk category
- Length: 240mm
- EN 455 parts 1, 2 & 3, AQL 0.65
- ISO 16604b:2003 viral test passed

**Supplied in packs of 100**

Code	Alt Ref	Size	Price
SAF5620	SL71	X small	£29.40
SAF5622	SL72	Small	£29.40
SAF5624	SL73	Medium	£29.40
SAF5626	SL74	Large	£29.40
SAF5628	SL75	X large	£29.40



## duoSHIELD™ PFT Latex 240 Gloves



- Natural rubber latex
- For general use in the laboratory and for patient care (dual risk)
- Powder free and fully textured with beaded cuff
- Ambidextrous
- Length: 240mm
- Class 1 Medical Device (MDR)
- Category III
- AQL 1.5
- EN ISO 374-1:2016 + A1:2018 Type B KPT (chemical risks)
- EN 16523-1:2015 + A1:2018 (chemical permeation)
- ISO 374-4:2019 (resistance to degradation by chemicals)
- EN ISO 374-5:2016 VIRUS (protection from viruses, bacteria and fungi)
- ISO 16604:2004 Procedure B (viral penetration test)
- EN 455-1/2/3/4:2000 (medical gloves for single use)

**Supplied in packs of 100**

Code	Alt Ref	Size	Price
SAF1048	654121	X small (6)	£39.70
SAF1050	654122	Small (7)	£39.70
SAF1052	654123	Medium (8)	£39.70
SAF1054	654124	Large (9)	£39.70
SAF1056	654125	X large (10)	£39.70