



Kevlar® Gloves



- 100% Kevlar® fibre gloves
- Resistant to cuts, abrasions, heat and flames
- Neptune coated – moisture resistant

Supplied in pairs

Code	Size	Price
SAF4040	Small	£18.73
SAF4048	Medium	£21.30
SAF4044	Large	£21.00



SENSITOUCH® Silk Glove Liners

- Effectively absorb perspiration helping to keep the hands dry and fresh
- Elimination of high levels of perspiration reduces the super hydration effect on the skin cells which causes them to swell and disrupt
- Positive barrier between the skin and the glove, helping to reduce any leaching of glove components which may cause skin irritation
- Ultra fine knitted silk fabric stretches to conform closely to the skin, resulting in virtually no loss of touch
- Absorption of perspiration increases operator tactility and work performance

Supplied in pairs

Code	Size	Cuffed	Price
SAF4440	Small (6.5 - 7)	No	£23.80
SAF4442	Medium (8 - 8.5)	Yes	£28.80
SAF4444	Large (8.5 - 9)	No	£23.80



Frosters Cryogenic Gloves

- For use with materials stored with liquid nitrogen (ranging from -120 to -190°C) and in other very cold environments
- Cut so that fingers and thumbs retain a high degree of movement
- Outer layer has a waterproof coating yet remains permeable to water vapour to prevent condensation inside
- Superb insulation needs less thickness giving more dexterity
- Inner layer of breathable cotton for comfortable wear
- Hanging loop inside the cuff
- Mid arm length: 380mm
- Category II
- EN 388:2016 2X2X (mechanical risks)
- EN 511:2006 22X (protection against cold)

Supplied in pairs

Code	Size	Price
CRY9000	Small (8)	£162.00
CRY9002	Medium (9)	£162.00
CRY9004	Large (10)	£162.00