

Tempshield® Cryo-Gloves®

- For ultra cold applications down to -196°C
- Maximum level of thermal protection, flexibility, and dexterity
- High performance
- Multi layered construction
- Thermal inner lining wicks moisture away from hands
- Category III
- EN 511:2016 (cold protection)
- EN 388:2016 (mechanical risks)

Supplied in pairs



Code	Length	Length, mm	Size	Price
SAF4060	Wrist	280 - 330 (wrist)	Small (8)	£279.00
SAF4062	Wrist	280 - 330 (wrist)	Medium (9)	£279.00
SAF4064	Wrist	280 - 330 (wrist)	Large (10)	£279.00
SAF4066	Mid arm	335 - 395 (mid arm)	Small (8)	£306.00
SAF4068	Mid arm	335 - 395 (mid arm)	Medium (9)	£306.00
SAF4070	Mid arm	335 - 395 (mid arm)	Large (10)	£306.00
SAF4072	Elbow	440 - 500 (elbow)	Medium (9)	£378.00
SAF4074	Elbow	440 - 500 (elbow)	Large (10)	£378.00
SAF4076	Elbow	440 - 500 (elbow)	X large (11)	£378.00

Heathrow Scientific Wire Glove Box Holders



- Keep gloves organised and ready for immediate use
- Free up benchtop space by mounting holders on the wall
- Sturdy, epoxy coated steel wire
- Fit most major brands of boxes
- Mounting screws included

Code	Capacity	Dims, w x d x h, mm	Price
SAF2150	1 glove box	140 x 108 x 209	£29.80
SAF2152	3 glove boxes	250 x 108 x 454	£61.00



View our Safety & PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

scientificlabs.co.uk/literature (UK & Northern Ireland)
scientificlabs.ie/literature (Republic of Ireland)
scientificlabs.com/literature (International)

