

3M™ 4440 Lab Coats Continued

3M Science.
Applied to Life.™

Code	Alt Ref	Size	Price
SAF0259	4440ZS	Small	£360.00
SAF0258	4440ZM	Medium	£360.00
SAF0257	4440ZL	Large	£360.00
SAF0260	4440ZX	X large	£360.00
SAF0254	4440Z2	2X large	£360.00
SAF0255	4440Z3	3X large	£360.00
SAF0256	4440Z4	4X large	£360.00



Kimtech™ A7 Lab Coats

KIMTECH

- Non-woven SMS, a secure and durable fabric with a cloth like feel
- Chemical and abrasion resistant spunbond 100% polypropylene outer layer, and a microfiber, cloth like inner layer designed to filter out a variety of water based liquids and dry particulates
- Triple stitched and overlapped seams
- Silicon and latex free
- Anti static
- Low lint
- Chemical and abrasion resistant spunbond 100% polypropylene outer layer, and a microfiber, cloth like inner layer designed to filter out a variety of water based liquids and dry particulates
- Category III
- EN 13034:2005 + A1:2009 Type 6 (limited chemical splash protection)
- EN 1149-1 (electrostatic properties)

Code	Alt Ref	Size	Price
SAF0030	9670	Small	£10.90
SAF0032	9671	Medium	£10.90
SAF0034	9672	Large	£10.90
SAF0035	9673	X large	£10.90
SAF0036	9674	2X large	£10.90



Kimtech™ A7 P+ Lab Coats

KIMTECH

- Non woven SMS fabric with triple stitched seams and abrasion resistant, film coated polypropylene
- Non slip outer layer and a cloth like feel
- Elasticated cuffs
- Thumb loops
- Extra length arms and high mandarin collars to enhance all round coverage and security
- Latex and silicone free
- Low lint
- High levels of particle and bacterial protection (with filtration efficiency of 99%)
- Individually packed and double bagged, making them suitable for ISO 7 (Grade C) environments
- Category III
- Particle shedding test (Helmke Drum) IEST-RP-CC003.3 category II
- EN 14126:2003 (infective agents)
- EN 1149-1 (electrostatic properties)
- EN 13034:2005 + A1:2009 Type 6 (PB) (limited chemical splash protection)

Supplied in packs of 15

Code	Alt Ref	Size	Price
SAF4768	97700	Small	£138.00
SAF4770	97710	Medium	£138.00
SAF4772	97720	Large	£138.00
SAF4774	97730	X large	£138.00
SAF4776	97740	2X large	£138.00