

Asecos CX Chemical Recirculating Filtration Cabinet for Non Flammable and Non Corrosive Chemicals Continued



Code	Alt Ref	Description	Dims, w x d x h, mm	Price
SAF6002	32933-050-31947	CX 4 door chemical cabinet with integrated recirculating extraction with 8 pull out shelves and sump tray	1055 x 520 x 2290	£6748.00
SAF6004	32952	Replacement filter - activated carbon	—	£714.00



Asecos FX Pegasus 90 Minute Fire Resistant Cabinet for Flammables



The FX range offers superior storage and filtration of flammables stored in accordance with BS EN 14470 part 1. The chemicals should be segregated in accordance with HSG71 and the safety data sheet. This unit can be used for flammable liquids and solids correctly segregated.

- Ninety minute fire resistance
- Choice of 3 widths: 1200, 900 and 600mm, all depths 615mm and all heights 2294mm
- Choice of standard door closing, automatic door closing
- Ready for connection to an extraction system
- Self closing doors, self damping and self sealing in the event of fire
- Choice of 7 door colours, anthracite grey surround
- Can be easily relocated within a building
- Packages with fixed shelves or drawers are available in either powder coated sheet steel or stainless steel
- Options include bottle racks, drawer mats and internal storage boxes
- Multi stage wideband filter, touch display monitor for service and user operation settings, air and fan monitoring, filter change required annually
- Ideal for where no on site extraction is available
- Delivery service – to goods in or point of use where the cabinet is delivered to where it is required, it is levelled and the packaging removed [site survey required for point of use]

Please enquire for more options

Code	Alt Ref	Description	Dims, w x d x h, mm	Standards	Price
SAF6006	32905-046-33635	FX Pegasus cabinet 120cm with extraction unit, 3 shelves, bottom sump and perforated insert	1200 x 615 x 2294	BS EN 14470	£7977.00
SAF6008	32950	Replacement filter - activated carbon	—	—	£811.00



Asecos G Ultimate 90 Minute Fire Resistant Cabinet for Compressed Gases



The G Ultimate 90 range offers the best storage for compressed gases stored indoors, in accordance with BS EN 14470 part 2. The compressed gases should be segregated in accordance with HSG71, British Compressed Gas Association guidance and the safety data sheet.

- 90 minute fire resistance
- Choice of 4 widths: 1400, 1200, 900 and 600mm, all depths 615mm and all heights 2050mm
- Includes rolling ramp and cylinders restraints
- Ready for mandatory connection to an extraction system
- Internal cabinet temperature kept to a maximum of 55 degrees in the event of fire
- Choice of 2 cabinet colours, warning yellow and light grey
- Bespoke options available for smaller cylinders
- Owing to the considerable weight of this range it is advisable point of use delivery, the cabinet is delivered to where it is required, it is levelled and the packaging removed [site survey will be required]

Please enquire for more options