

## Asecos G Ultimate 90 Minute Fire Resistant Cabinet for Compressed Gases Continued

Code	Alt Ref	Description	Dims, w x d x h, mm	Standards	Price
<b>SAF6010</b>	30678-001-30682	G Ultimate 90 60cm width cabinet with rolling ramp, cylinder restraints and easy connection to pipe lead throughs, cabinet can hold 2 x 50 litre cylinders	600 x 615 x 2050	BS EN 14470 part 2	<b>£4208.00</b>

## Asecos Under Bench 90 Minute Fire Resistant Cabinet for Oxidising Substances

The UB range offers a broad range of under bench storage for flammable liquids, it is manufactured to BS EN 14470 part 1, the cabinets are ideal for storage where space is premium. They have proven themselves in many modern laboratory designs, whether building a new laboratory or refurbishing an existing laboratory they are the perfect solution. Contents must be stored in line with HSG71 and the chemical safety data sheets.

- Ninety minute fire resistance
- Choice of widths: 1400, 1100, 890 and 593mm, choice of depths 500 and 574mm, standard height 600mm, 800mm height available only in 1100 and 593mm widths, also with wing doors or pull out draw front
- Ready for connection to an extraction system through a 50mm port
- Self closing doors, self damping and self sealing in the event of fire
- Light grey cabinets
- Internal options of second level draws and fixed shelves
- All cabinets contain a sump and perforated insert avoiding displacement of sump volume in the event of a spill
- Delivery service – to goods in or point of use where the cabinet is delivered to where it is required, it is levelled and the packaging removed



Please enquire for more options

Code	Alt Ref	Description	Dims, w x d x h, mm	Standards	Price
<b>SAF6020</b>	30472-001-35286	Under bench cabinet with pull out draw front and perforated insert	593 x 600 x 600	BS EN 14470 part 1	<b>£2004.00</b>

## Asecos Q Classic 30 Minute Fire Resistant Cabinet

The Q Classic range offers superior storage of flammable liquids stored in accordance with BS EN 14470 part 1. The liquids should be segregated in accordance with HSG71 and the safety data sheet.

- 30 minute fire resistance
- Choice of 3 widths: 1164, 864 and 564mm, all depths 560mm and all heights 1947mm
- Choice of standard door closing, automatic door closing
- Ready for connection to an extraction system
- Self closing doors, self damping and self sealing in the event of fire
- Choice of 7 door colours, anthracite grey surround
- Can be easily relocated within a building
- Packages with fixed shelves or drawers are available in either powder coated sheet steel or stainless steel
- Options include bottle racks, drawer mats and internal storage boxes
- Delivery service – to goods in or point of use where the cabinet is delivered to where it is required, it is levelled and the packaging removed (site survey required for point of use)



Please enquire for more options